




TV & Audio



Bluetooth Audio



Conversation



Home Alerting



Alarm Clocks



Telephones



Tinnitus



Commercial Audio



www.sarabec.com

About this Catalogue

SARABEC® AIMS TO MAKE LIFE A LITTLE EASIER FOR THOSE WITH HEARING DIFFICULTIES BY PROVIDING HIGH QUALITY PRODUCTS THAT HELP YOU:

- Stay in touch with friends and relatives
- Know what's happening around the house
- Enjoy recreational time without feeling isolated

We have over 40 years of experience in providing products that offer simple hearing solutions for use both with or without a hearing aid.

We aim to supply all of our customers, both private individuals and companies alike, with a high quality service that's appropriate to their needs.

WE ARE PLEASED TO OFFER YOU OUR LATEST SELECTION OF PRODUCTS TO HELP YOU WITH:


- Listening - enjoy TV, radio, and conversation at a volume suitable for you and others
- Alerting - doorbells, telephones, and alarms
- Telephones - amplified phones including corded, cordless, and mobile

BUYING FROM SARABEC

- Easy ordering online or over the phone
- Technical help when you need it - both before and after sales

1 MONTH TRIAL

With our 1-month trial policy, you can be sure that the products you choose are right for you by trying them in



the comfort of your own home. If you are not happy with the product, simply contact us and we will let you know the best way to return the product to obtain a refund. We cannot give you a refund if goods, when returned, are found to be damaged, incomplete or have been tampered with.

TECHNICAL HELP

If, at any time, you wish to discuss whether a product is right for you before ordering, then contact us, and our friendly, knowledgeable staff will be happy to help - they can also support you after your purchase if you need assistance when using your product.

VAT RELIEF

If you have hearing loss and are buying for your own use, you can obtain VAT relief on certain products in this catalogue. Products that are eligible for VAT relief are indicated as shown below.



This product does not qualify for VAT Relief.



VAT Relief is available on this product.

To obtain VAT Relief you must complete a VAT Relief Certificate at the time of your order. When ordering online simply tick the relevant section. VAT refunds cannot be given afterwards. If you are unsure, or want to know more, please contact us.

Note: Some accessories, if not bought as part of a system, cannot be given VAT Relief.

EVERYDAY LIFE WITH A PACEMAKER

Pacemaker manufacturers give detailed instructions about possible sources of interference from products that you should avoid. Please refer to your booklet about your type of pacemaker. If you have any further questions about what you can and can't do, please contact your pacemaker clinic. Further information is given on our website.

SYMBOLS USED THROUGHOUT THIS CATALOGUE:

| | |
|--|---|
| A blue circular icon containing a white 'T' symbol. | This product is suitable for users with hearing aids that have the "T" Loop setting or Loop Program |
| A yellow circular icon with a black border containing the text "WARRANTY 2 YEAR WARRANTY". | This product comes with a warranty for the number of years shown |
| A blue circular icon containing the white Bluetooth symbol. | This product uses Bluetooth connectivity technology |

SARABEC LIMITED,
15 HIGH FORCE ROAD,
MIDDLESBROUGH,
TS2 1RH, UK

OFFICE HOURS: 8.15am to 4.45pm - Mon-Thur | 8.15am to 4.00pm - Friday
TELEPHONE: 01642 247789 | **FAX:** 01642 230827 | **TEXTPHONE (NGT):** 18001 07497 242 294
E-MAIL: mail@sarabec.co.uk | **WEBSITE:** www.sarabec.com

Contents

03 TV & AUDIO

- 03 Digital Wireless TV Listening Systems
- 08 Induction Loop Systems
- 10 TV/Audio/Loop Accessories

Have difficulty hearing your television, radio, Hi-Fi or other audio? Bring the wanted sound direct to your ears without disturbing others. Our systems are designed to be used with or without hearing aids. Explore our portable & stationary systems, & enjoy recreational time again.

13 BLUETOOTH AUDIO

- 14 Equipment For Mobile Phones
- 15 Personal & TV Listening
- 16 Sleep Headphones & Sound Players

More and more hearing aids have Bluetooth capability built into them. You can now connect hearing aids to other bluetooth devices such as phones or TV transmitters via the "T" Loop setting or Loop Program.

17 PERSONAL LISTENING & COMMUNICATION

- 18 Listening Options
- 19 Personal Listeners
- 23 Personal Audio Accessories

Designed to be used with or without hearing aids, our Listening & Conversation products aim to make life easier if you have difficulty in daily listening situations or with conversations at home or otherwise. Versatile and portable, these systems use a range of listening options to bring sound closer to you.

25 HOME ALERTING

- 26 Signolux Home Alerting
- 35 Fireangel Smoke/CO Alarms
- 37 Doorbells & Accessories
- 42 Geomarc Amplicall Baby Monitor
- 43 Alarm Clocks

These products make you aware of everyday occurrences in the home or at work via vibrations and flashing lights, so you can "feel" and "see" the alerts. Get alerted to the ringing of the telephone or doorbell, the sound of the smoke or fire alarm or a baby crying.

47 PHONES & ACCESSORIES

- 49 Mobile Phones
- 51 Portable DECT Telephones
- 58 Corded Desktop Telephones
- 63 Phone Accessories

If telephone callers' voices seem very quiet or muffled, you will find a solution amongst our selection of telephones and amplifiers. Explore page 63 to view products that let you know the phone is ringing.

64 TINNITUS

- 65 Wellcare & Sound Oasis

Manage Tinnitus symptoms and reduce the effects of buzzing/ringing in your ears with our tinnitus relief products that create a relaxing atmosphere, reduce stress, and help you sleep.

67 COMMERCIAL AUDIO

- 68 LA90 Portable Induction Loop
- 69 ProLoop Systems
- 70 Accessories

We can work with you to improve communication for hearing impaired people throughout non-domestic environments like in a church, village hall or any public space.

71 ORDER FORM & INFORMATION

You can complete your order via this form and attach the VAT Relief Certificate (if applicable), and send it back to us in the post.



TV & AUDIO

DIGITAL & WIRELESS TV LISTENING SYSTEMS

WHAT IS A TV LISTENER?

TV listeners let you listen to your audio devices (TV, radio or other device), at a volume that you are comfortable with, without disturbing others. With no trailing wires, TV listeners bring the sound directly to your ears.

IS A TV LISTENER RIGHT FOR ME?

Do you want to listen to the TV at a higher volume than other members of the family? Is listening in a noisy situation stressful? Is turning the volume up on the TV or radio not enough? Do you have difficulty in hearing the programmes clearly?

If friends or family complain that your TV or radio is too loud then that's when one of our TV Listening systems may prove invaluable as they can drastically improve your entertainment enjoyment. With these wireless TV listening systems, you can choose how loud you want the sound to be without disturbing others. You can also move around freely without the need to have trailing wires between you and the device you're listening to (Note: you will still need to connect the

transmitter to your TV, radio or music to pick up the sound using the leads provided).

TV LISTENERS FOR HEARING AID USERS

Your hearing aid may have the "T" Loop setting or Loop Program. Select one of the products with the T symbol next to it to make the most out of your hearing aid and enhance your listening. Sounds will be transmitted directly to your hearing aid's T loop program, increasing intelligibility and significantly reducing background noise. Products without the T symbol are suitable for those who do not wear hearing aids. All TV Listeners require you to wear a receiver whether or not you use a hearing aid.

HOW DO TV LISTENERS WORK?

Digital listening systems use the 2.4GHz frequency to transmit their sound from the transmitter to the receiver. They use a coding system that means they are unlikely to suffer any form of interference. The coding system requires the transmitter to be matched/paired to the receiver(s) required. Multiple receivers can normally be used with a single transmitter.



sonumaxx^{NX}

SONUMAXX - A SIMPLE-TO-USE TV LISTENING SYSTEM FOR BOTH HEARING AID AND NON-HEARING AID USERS.

With its clear sound, the latest Sonumaxx Digital Wireless System will help you listen more clearly to your TV, radio, or other sound sources without any trailing leads.

With the Sonumaxx system you can switch between the transmitter and internal microphone in the receiver when listening to the TV sound. This is because Sonumaxx receivers have a built-in communication microphone and amplifier function, this allows the user to receive amplified sounds from conversations all at the press of a button on the receiver, very useful one to one speech with friends or family!

When listening to the TV, you can listen at your own volume without disturbing others. You can also move around whilst listening to TV, radio or Hi-Fi as the signal is not interrupted by walls. The system has a range of up to 25 metres and the headset systems have an output of up to 125dB.

OUR RANGE OF SONUMAXX SYSTEMS

There are two types of systems available; the standard Headset System where the headset and receiver unit are a single device, and the PR systems; these allow you to plug any audio listening accessory into the receiver unit such as headphones, neckloop, ear hooks etc. You can see all 4 options here.

SYSTEM COMPONENTS

All systems consist of two parts; a mains powered transmitter and a receiver powered by a rechargeable battery. The transmitter picks up the sound from your TV using the supplied connection lead kit or microphone. This sound is then sent to the receiver. The Sonumaxx systems are already paired for your convenience and to avoid interference.

BATTERY LIFE

A fully charged battery can give up to eight hours listening time depending on the volume level used. The system is supplied with a second battery to provide continuous listening. An intelligent charging system ensures batteries will be recharged fully in less than 3 hours.

Sonumaxx NX Headset System

Transmitter, Conventional Under-chin Headset Receiver & Connection Lead Kit with Microphone.



Sonumaxx NX Headset System

Order code **230-416**
Price **£119.90 Ex. Vat**
£143.88 Inc. Vat

Colour **Black/Grey**
Size **Base: 125mm Ø x 40mm**
Receiver: 225x125x 25mm



Sonumaxx NX PR Headset System

Transmitter, PR Pocket Receiver, Plug-in Stethoscope Headset & Connection Lead Kit with Microphone.



Sonumaxx NX PR Headset System

Order code **230-425**
Price **£119.90 Ex. Vat**
£143.88 Inc. Vat

Colour **Black/Grey**
Size **Base: 125mm Ø x 40mm**
Receiver: 225 x125x25mm



Sonumaxx NX PR Neckloop System

Transmitter, Pocket Receiver, Plug-in Neckloop & Connection Lead Kit with Microphone.



Sonumaxx NX PR Neckloop System

Order code **230-426**
Price **£119.90 Ex. Vat**
£143.88 Inc. Vat

Colour **Black/Grey**
Size **Base: 125mm Ø x 40mm**
Receiver: 225 x125x25mm



Sonumaxx NX PR Headset/Neckloop System

Transmitter, Pocket Receiver, Headset, Neckloop & Connection Lead Kit with Microphone.



Sonumaxx NX PR Headset/Neckloop System

Order code **230-427**
Price **£131.90 Ex. Vat**
£158.28 Inc. Vat

Colour **Black/Grey**
Size **Base: 125mm Ø x 40mm**
Receiver: 225 x125x25mm



TV & Audio

earis®

TV & Audio Listening

Telephone (voice): 01642 247789
Fax: 01642 230827 Textphone (NGTS): 18001 07497 242 294

www.sarabec.com
Email: mail@sarabec.co.uk

05

SARABEC® - Hearing products for a better lifestyle

PREMIUM
Digital Wireless Listening



EARIS PREMIUM WIRELESS SYSTEMS

WHAT IS EARIS PREMIUM WIRELESS LISTENING?

The Earis System is a premium, digital TV listening system that uses dynamic sound technology to deliver the highest sound quality and clarity that is currently available in TV listening equipment today. Brilliant for people with or without hearing aids to watch TV or just listen to music. The system is easy to use and is ready to connect to TV and audio systems.

HOW DOES THE EARIS PREMIUM WORK?

Earis significantly improves understanding of speech in audio even if background sounds are present.

You can choose from 5 listening profiles for a personalised listening experience.

The fast-acting 'Automatic Sound Profiling' feature eradicates sudden volume changes in sound after you've selected your chosen listening level - great for when changing channels or when adverts cut in.

EXCELLENT 1:1 COMMUNICATION

Earis receivers have a communication amplifier function, where the built-in microphone in the Earis receiver amplifies sound in the same room, making it excellent for 1:1 communication. Earis then cleverly stores this setting to ensure that the user always returns to the last set volume when switching over.

SELECT FROM TWO EARIS OPTIONS:

Two Earis Premium systems are available; a headset

receiver system, and a PR Neckloop System. Each system uses the same transmitter, which has connections and cables for digital audio and analogue audio signals, both of which are usually present on all TV's and also Sky, Virgin and Freeview boxes. The PR Neckloop System allows you to plug in your choice of listening device (headset/earphones/ear hooks/or anything with a standard 3.5mm jack plug).

THE SYSTEMS FEATURE:

- Built-in receiver charging cradle and charging port for an optional second battery
- Connection kit for TV and audio systems including an optical digital audio lead, audio connection lead with 3.5mm stereo plugs, and adapter to RCA phono plugs
- Up to 122dB amplification
- Replaceable Li-Pol rechargeable battery, with fast full charge in under 2 hours - 80% in 30 minutes
- Approximately 8 hours listening time per full battery charge
- Large tactile volume wheel with LED level display
- Easy-to-use communication amplifier available at the press of a button on the receiver
- Tone & Balance controls
- PR Pocket Receiver: Output connector 3.5mm stereo jack socket, 16ohm minimum

Through extensive research into the listening needs of hearing impaired people, it has been determined that Earis's 5 listening profiles match the needs of 95% of users when they are listening to speech and/or music from film, TV, radio or other sound source. Earis delivers constant high-quality sound using the selected profile.

Earis Premium Digital Headset System

Earis Premium Digital Headset System Transmitter, Under-chin Headset Receiver & Connection kit. Amplification level (up to 122dB+).



Earis Premium Digital Neckloop System

Earis Premium Digital PR Neckloop System Transmitter, Receiver, Plug-in Neckloop & Connection kit.



Earis Premium Digital Headset System

| | | | |
|------------|---|--------|----------------------------------|
| Order code | 230-190 | Colour | Graphite/Silver |
| Price | £249.00 Ex. Vat £298.80 Inc. Vat | Size | Headset: 244x112x24mm |



Earis Premium Digital Neckloop System

| | | | |
|------------|---|--------|---------------------------------|
| Order code | 230-191 | Colour | Graphite/Silver |
| Price | £264.00 Ex. Vat £316.80 Inc. Vat | Size | Receiver: 63mmØx24mm |





CL7500 OPTI WIRELESS HEADPHONES




THE GEEMARC CL7500 OPTI IS DESIGNED FOR HEARING-IMPAIRED PEOPLE.

It has an adjustable volume of up to 125dB SPL, making it suitable for people with hearing loss ranging from low to severe. The headset is lightweight, weighing in at only 260 grams. The earpieces contain controls which include Left & Right Balance, & volume and tone. You can charge your wireless TV headset by placing it on its base, using an external battery, or even connecting it directly to your laptop! It takes 6 hours for the batteries to fully charge.

The CL7500 main base can be plugged into your televisions digital optical or analogue output sockets using the leads provided. It also works with other audio devices such as radios & DPD players, as long as those devices have a compatible audio output socket.

CL7500 OPTI Wireless Headphones

| | | | |
|------------|-------------------------------------|--------|---|
| Order code | 252-221 | Colour | Silver/Black |
| Price | £141.66 Ex. Vat £169.99 Inc. Vat | Size | Base: 205x130x75mm Headset: 190x160x65mm |



- ✓ Built-in digital treble/bass adjustment (+/-10dB)
- ✓ Amplification up to +125dB
- ✓ Balance Control (left/right)
- ✓ Digital optical & analogue connections

Telephone (voice): 01642 247789
 Fax: 01642 230827 Textphone (NGTS): 18001 07497 242 294

www.sarabec.com
 Email: mail@sarabec.co.uk


WHAT IS AN INDUCTION LOOP SYSTEM?

Induction loop systems transmit sound directly to a hearing aid's loop ("T") facility or program which is built into or can be enabled in most hearing aids.

The main advantage of using the loop facility is that it can minimise the distracting effect of background noise and so enhance the sounds you want to hear making them easier to understand. Because the loop receiver facility is built into your hearing aid you will not need to wear an additional receiver.

A small amplifier is positioned near the television or sound source, it is then plugged into the television/sound source using one of the supplied direct connection cables (most commonly used is the optical digital cable - the socket for this may be marked up on your TV as Optical, SPDIF, TosLink or Digital Audio).

A very thin wire runs around the perimeter of the room to and from the amplifier. This wire can be positioned underneath floor surfaces such as carpet, laminate, or wood. It can also be attached to/follow your skirting boards or picture rail, or even be placed in a ceiling cavity or attic.



Some users find that though they benefit from the sound of the TV being directly transmitted to their hearing aid, they feel

cut off from other useful sounds like the doorbell or their partner speaking. This can be overcome by asking your audiologist for an additional program where the loop and microphone are mixed together (sometimes known as M/T or mixed loop and mic).

The real beauty of a loop system is that once installed, and your chosen volume and tone levels are set to how you like them, there is little to no interaction required with the amplifier.

Not sure how to use the loop features of your hearing aid? Ask your audiologist as you may need to activate the features, or adjust the sensitivity of the loop program.

Both the LA240 Room Loop System and the LA240 Loop Pad are supplied with lead kits that let you connect the amplifier directly to your TV with digital optical, phono and jack connectors, and they both come with a microphone that has a 1.2m lead (can be used for conversation or monitoring the telephone or doorbell ring). See next page for other connection possibilities.

DIFFERENCES BETWEEN THE ROOM LOOP SYSTEM AND LOOP PAD SYSTEM:

The room loop system will cover an entire room, and you will pick up the sound from the system wherever you are within the room. However a loop pad system generates a very localised loop signal via its small loop pad which goes under a cushion or on the back of a chair. You will only pick up the signal while sitting in that one place.

LA240 Room Loop System

This system will cover an entire room.



LA240 Room Loop System

40 metres (approx 130ft) of thin, white loop cable (2.5mm diameter) and cable clips

Order code **135-250**
Price **£134.90 Ex. Vat**
£161.88 Inc. Vat

Colour **Silver**
Size **180x140x42**



LA240 Loop Pad System

Ideal for a single person or where a limited loop field is required. (Supplied with 5m (15ft) connecting lead.



LA240 Loop Pad System

This LA240 system includes a pad in place of the cable reel.

Order code **135-255**
Price **£134.90 Ex. Vat**
£161.88 Inc. Vat

Colour **Silver**
Size **180x140x42**



NEED HELP CONNECTING OUR EQUIPMENT TO A TELEVISION OR OTHER AUDIO SOURCE?

CONNECTING DIRECTLY TO A TELEVISION

When connecting directly to a TV set via optical, you must ensure the TV is set to PCM. To do this, go into the audio/sound settings menu on the TV and look for 'Digital Audio Output/Optical Out', normally this would be set to 'Dolby Digital', it needs to be changed to 'PCM'. Please note that some televisions auto switch this setting!

NOTE: Using the microphone supplied with your equipment attached to the speaker grille of the television will always return good results and is an excellent back-up way to connect to the TV.

SETTING UP YOUR EQUIPMENT WITH DIRECT CONNECTION TO SKY Q & SKY HD BOXES VIA OPTICAL

- » Press **Home** on your Sky Q remote, select **Settings**, then select **Setup**.
- » Select **Audio visual**, then select **Digital Audio Output** optical and set it to **NORMAL**
- » Press the **GREEN** button to save the settings
- » Power **OFF** the Sky Box, leave for 30 seconds and turn back on, leave for 5 minutes so the box can re-calibrate itself and then test the system.

Note: Sky Glass TV's & Sky Stream Pucks are not compatible with our equipment via a direct connection. However, a microphone may still be used. Please contact us for further details if required.

SETTING UP EQUIPMENT WITH DIRECT CONNECTION TO VIRGIN MEDIA SET TOP BOXES VIA OPTICAL





- » Press **Home** on your remote
- » Select **Help & Settings**
- » Select **Settings**
- » Select **Audio**
- » Choose **Dolby Digital to PCM**. This setting applies to both HDMI and optical out.
- » Power down the box for 5 minutes, restart the box, leave for a few minutes and test the system
- » **Virgin Media 360 Boxes ONLY: Select Settings > Audio & Video > PCM2.0> Restart Box as above.**

WHY DON'T I GET AUDIO FROM DIGITAL SOURCES (NETFLIX, APPLE TV, FIRE STICK, ETC)?

Audio coming from certain content providers or devices (Netflix, Amazon Prime, Apple TV etc.) is set to output Dolby and not PCM by default. To resolve this, you will need to change the digital audio output format in the app or device to PCM, firstly try the normal audio settings within the app. Usually, these audio settings are only available once something is playing on screen. This does NOT guarantee that you will hear the audio from these sources through your device but it does give the best chance.

N.B When connecting directly to an external device such as SKY or Virgin Media you will not hear channels and smart services from the TV through the connected equipment.

COMMON SOCKET TYPES

| | |
|--|---|
| <p>Optical/SPDIF/Digital Audio/TOSLink Preferred connection method, will generally require a change to the TV or connected equipment to PCM audio output.</p> |  |
| <p>HDMI Can be used in conjunction with the Sarabec ARC Audio Converter.</p> |  |
| <p>PHONO/RCA Can be a good option, but make sure that the RCA sockets are marked as Audio Output and NOT Audio Input.</p> |  |
| <p>3.5mm Audio/Headphone A good method of connecting. However, in most cases, using this socket will turn off the internal speakers of the television.</p> |  |



DIGITAL CONNECTION CONVERTERS

Dolby Digital Audio Decoder Kit

A clever little device which will, in most cases, allow you to connect any of our Television Listening Systems (LA240, Sonumaxx & Earis) to any television, set top box, or audio output device without needing to change the settings on the television to output PCM signals. This is especially useful when the television or other audio source doesn't auto switch to PCM output or when you cannot find the appropriate settings in the menu systems.



| Dolby Digital Audio Decoder Kit | |
|---------------------------------|------------------------|
| Order code | 354-132 |
| Price | £29.99 Inc. Vat |
| Colour | Black |
| Size | 105 x 75 x 25mm |



Sarabec HDMI ARC Audio Extractor Kit

Connect your existing listening equipment to your television's ARC HDMI socket with the Sarabec ARC Audio Extractor Kit. Many televisions do not have an optical or traditional audio output sockets, so using this adapter will allow you to extract audio from the TV. This is especially useful if the television doesn't have optical or traditional audio sockets.

It can also be used in place of a traditional digital to analogue converter. Supplied with 0.5m HDMI lead and USB power cable.



| HDMI ARC Audio Extractor Kit | |
|------------------------------|------------------------|
| Order code | 354-130 |
| Price | £39.98 Inc. Vat |
| Colour | Black |
| Size | 100 x 70 x 25mm |



AUDIO & VIDEO MISC.

Loop Pad 340 X 340mm

Use with loop amplifiers and systems if the area to be looped needs to be small and discreet. Use at home, place on a chair or chair back, and in reception areas under counters. Also used to demonstrate loop systems to clients.



Loop Pad 340 X 340mm

| | |
|------------|---|
| Order code | 110-053 |
| Price | £24.99 Ex. Vat £29.99 Inc. Vat |

| | |
|--------|------------------|
| Colour | Black |
| Size | 340x340mm |



Standard Microphone 1.2m lead

Standard microphone for many of our products including amplifiers.



Standard Microphone 1.2m lead

| | |
|------------|------------------------|
| Order code | 330-005 |
| Price | £17.99 Inc. Vat |

| | |
|--------|--------------|
| Colour | Black |
|--------|--------------|



40m Single Core Loop Cable

A 40m reel of single-core loop cable for use with our current LA240 loop amplifier as well as older LA100/LA200/LA110/LA210 & LA215 home loop amplifiers.



40m Single Core Loop Cable

| | |
|------------|------------------------|
| Order code | 340-037 |
| Price | £17.88 Inc. Vat |

| | |
|--------|---|
| Colour | White |
| Size | Length - 40 metres 0.75mm2 diameter Impedance - 1.5 ohms |





ACCESSORIES



Geemarc Easy-to-use TV1 Remote Control

The TV1 remote control is very easy to use with its simple and large design. The channel and the volume can easily be adjusted with the large programmable buttons.

- Extra-large buttons for easy use
- Equipped with 19 programmable buttons

Geemarc Easy-to-use TV1 Remote Control

| | | | |
|------------|---|--------|--------------------|
| Order code | 252-222 | Colour | Black |
| Price | £10.83 Ex. Vat £13.00 Inc. Vat | Size | 166x46x25mm |



Doro HandleEasy 321rc Remote Control



The HandleEasy 321rc is specially designed for people who have limited dexterity. It is very easy to hold and to use for changing channels and adjusting the volume of your TV or home stereo.

- 7 programmable keys
- Hard-keys for best tactile feeling
- Learning function - learn commands from your existing remote controls (program channel keys for a set-top box and volume control for a TV)
- Compatibility to RC5/RC6 IR equipment (all major brands of TV/Video/DVD/stereo/set-top boxes). Includes dual-command support
- 3 x AAA batteries (included)



Doro HandleEasy 321rc Remote Control

| | | | |
|------------|------------------------|--------|--------------------|
| Order code | 253-140 | Colour | White/Grey |
| Price | £24.98 Inc. Vat | Size | 165x52x30mm |



Bluetooth Audio



WHAT IS BLUETOOTH?

Bluetooth is an easy to use short-range wireless communications technology now found worldwide in millions of products that we use every day - including headsets, smartphones, laptops, and portable speakers. Bluetooth technology is used for short range wireless communication and listening. It is a secure system that allows you to enjoy sound that is free of interference and the distractions of background noise. Bluetooth devices radiate much less power than a mobile phone and typically have a useful range for most devices between 5 and 10 metres.

BLUETOOTH TECHNOLOGY FOR HEARING AIDS

Bluetooth is increasingly being used to connect hearing aids to other devices. Hearing aid streamers often use Bluetooth to connect to phones or TV Bluetooth transmitters and then send those signals to the hearing aid via the aids own radio system or "T" Loop setting or Loop Program. More and more hearing aids have Bluetooth capability built into them.

BLUETOOTH NECKLOOPS

Bluetooth neckloops are universal receivers, and unlike streamers, can be used with any hearing aid that has a "T" Loop setting or Loop Program. So, if you change your brand of hearing aid, you won't need a new streamer. They can be paired with mobile phones and any other device that has a bluetooth output. Bluetooth neckloops allow you to freely move around whilst being connected to your mobile phone, even if it's sitting on your desk, left on the couch or in your bag. They enable discreet conversation without using the hands-free mode.

BLUETOOTH® ADVANCEMENTS

As technology improves, it's now possible to connect (pair) more than one device with the Bluetooth receiver, so switching between a phone and TV becomes easier. More recent introductions of Bluetooth (current version: 5.0), means that the delays in receiving sound, that some people previously experienced, have been rectified, including delays in experiencing lip-sync, and more importantly, the power required to run the receivers.

BLUETOOTH EQUIPMENT FOR MOBILE PHONES

KEEP YOUR HANDS FREE WHILST MAKING AND ANSWERING PHONE CALLS.

The CM-BT2 Bluetooth Neckloop uses short range Bluetooth wireless technology, giving you a wireless connection between your Bluetooth-enabled phone and your hearing aid via the "T" Loop setting or Loop Program.

Although primarily designed for use with mobile phones, there are plenty of other devices you can now communicate with when using a Bluetooth neckloop.



The sound is received by the CM-BT2 and transmitted direct to the loop receiver in your hearing aid. You can control the volume of the loop directly on the receiver. A built-in microphone allows full two-way communication.

- Power supply and charging cable
- 3.5mm socket to listen via other accessories

CM-BT2 Bluetooth Neckloop

Order code **230-007**
 Price **£84.99 Ex. Vat**
£101.99 Inc. Vat

Colour **Black/Blue**
 Size **48x76mm**



CM-BT2 BLUETOOTH NECKLOOP



HUMANTECHNIK



Artone 3 Max Bluetooth Loopset

Artone 3 MAX streamer is the most powerful bluetooth neckloop on the market, it works even for those with severe hearing loss and cochlear implants. You will not need Bluetooth hearing aids. This will connect your Bluetooth devices to a standard loop enabled hearing aid.

When using the Artone 3 Max, your loop-enabled hearing aid will have Bluetooth connectivity, this will allow the user to communicate using a mobile phone as a hands-free device. You will also be able to listen to music and the TV on other Bluetooth devices via hands-free and wire-free. Can also be used with other equipment in the Artone range.

Artone 3 Max Bluetooth Loopset

Order code **233-005**
 Price **£109.90 Ex. Vat**
£131.88 Inc. Vat

Colour **Black**
 Size **45x30x15mm**



Bluetooth Audio



BLUETOOTH EQUIPMENT FOR PERSONAL & TV LISTENING

PHONAK
life is on



Phonak TV Connector

The Phonak TV Connector is a wireless hearing aid accessory that streams TV and music directly to your hearing aids, transforming your Bluetooth-enabled hearing aids into wireless stereo headsets. A simple plug-and-play interface works with TVs and other audio sources to support direct hearing aid connectivity to Lumity, Paradise, and Marvel hearing aids.

| Phonak TV Connector | | | |
|---------------------|-------------------------|--------|-------------------|
| Order code | 240-127 | Colour | Black |
| Price | £159.00 Ex. Vat | Size | 50x50x50mm |
| | £190.80 Inc. Vat | | |

Artone 3 Max Bluetooth Loopset

Artone 3 MAX streamer is the most powerful Bluetooth neckloop on the market - it works even for those with severe hearing loss and cochlear implants.

You will not need Bluetooth hearing aids. The Artone 3 Max will connect your Bluetooth devices to a standard loop-enabled hearing aid, giving them Bluetooth connectivity, & providing you with a hands-free/wire-free experience when using your mobile phone, or when using other Bluetooth devices to listen to music and TV.



| Artone 3 Max Bluetooth Loopset | | | |
|--------------------------------|------------------------|--------|-------------------|
| Order code | 233-005 | Colour | Black |
| Price | £109.90 Ex Vat | Size | 45x30x15mm |
| | £131.88 Inc Vat | | |



Artone MIC Wireless Bluetooth Microphone

A wireless microphone that connects to your Artone Bluetooth Loopset. It will help you to hear the sounds your were missing until now.

The wireless microphone can be used by someone speaking at a distance from you and where there is background noise.

- Stereo frequency range - A2D2
- Microphone sensitivity control
- Direct Audio Input option
- Rechargeable
- Microphone with USB charging cable
- Tie clip microphone
- 10cm 3.5mm jack to jack lead



| Artone MIC Wireless Bluetooth Microphone | | | |
|---|-------------------------|--------|---------------------|
| Compatible with Artone 3 Max & many other Bluetooth receivers | | | |
| Order code | 233-020 | Colour | Black/Silver |
| Price | £85.40 Ex. Vat | Size | 50x15x10mm |
| | £102.48 Inc. Vat | | |

Artone TVB TV Transmitter (Optical)

The TVB transmitter only works in conjunction with the Artone 3 Max. If you can hear the TV but with all the background noise of your surroundings - use TVB with the Artone 3 Max. It will transfer the sound from your TV direct to your loop-enabled hearing aid without interference.

- Transmitter with power cable
- Optical Audio Lead
- 10cm 3.5mm jack to jack lead



| Artone TVB TV Transmitter (Optical) | | | |
|---------------------------------------|------------------------|--------|---------------------|
| Supplied with 2 x 3.5mm Phono Adapter | | | |
| Order code | 233-022 | Colour | Black/Silver |
| Price | £49.90 Ex. Vat | Size | 60x60x60mm |
| | £59.88 Inc. Vat | | |



BLUETOOTH SLEEP HEADPHONES & SOUND PLAYERS

SNOOZEBANDS EASE TINNITUS SYMPTOMS THROUGH BLUETOOTH CONNECTIVITY.

Compatible with **any** Bluetooth device, these Snoozebands connect to your mobile phone so you can listen to relaxing sounds to ease your tinnitus when trying to sleep or even just use it to listen to music or audio books. The Snoozebands range from 10-15 hour playtime, offer soft padded speakers and a built-in microphone to answer calls. Simply connect to your smartphone or Bluetooth player and use your favourite apps or sounds to help with your tinnitus. Soft, breathable, comfortable, and washable material. There are 3 Snoozebands to choose from.

SnoozeBand™ Original Grey

Enjoy a deeper, more comfortable sleep with the original SnoozeBand™. The ultra-soft thin padded speakers offer unmatched comfort and an immersive audio experience. It's perfect for side sleepers. Includes integrated microphone. One size fits all.



- Automatic shut-off after 1 or 2 hours
- Up to 10-hour playtime from 1 charge
- No disturbing flashing light

SnoozeBand™ Original Grey

| | | | | | |
|------------|------------------------|--------------|-------------------|--|--|
| Order code | 260-100 | Colour | Grey | | |
| Price | £38.99 Inc. Vat | Battery life | 8-10 hours | | |

SnoozeBand™ Pro 2.0

The updated SnoozeBand™ Pro model, with thinner, built-in padded speakers and a 3D eye mask for a true black-out experience. Made of soft brushed cotton on the outside and lightweight padding around the eyes, allowing you to relax & enjoy your favourite sounds.



- Blocks 100% of light
- 1 charge (up to 10 hours of playtime)
- Strap to adjust size

SnoozeBand™ Pro 2.0

| | | | | | |
|------------|------------------------|--------------|-------------------|--|--|
| Order code | 260-104 | Colour | Black | | |
| Price | £44.99 Inc. Vat | Battery life | 8-10 hours | | |

SnoozeBand™ Deluxe

The most advanced sleep headphones on the market. Ultra-thin (4mm) speakers covered in foam for unparalleled side-sleep comfort. Improved sound quality & longer battery life. Portable, blocks out 100% of light, has adjustable headphone position & strap, and features a built-in microphone for calls.



- Up to 15-hour playtime from 1 charge
- Built-in zip for ice & heat pack

SnoozeBand™ Deluxe


| | | | | | |
|------------|------------------------|--------------|--------------------|--|--|
| Order code | 260-107 | Colour | Grey | | |
| Price | £49.99 Inc. Vat | Battery life | 13-15 hours | | |

SOUND OASIS BLUETOOTH SPEAKER BST-100

Sleep soundly with any one of the 10 built-in sounds to help relieve tinnitus and to provide a calming background noise for relaxation.

Perfect for improving sleep, blocking noise, managing tinnitus and calming babies. Connect to a mobile phone or computer wirelessly and stream additional sounds, music or even an audio book. Alternatively, connect via an auxiliary cable (not supplied).

Exceptional big bass sound in super compact size, perfect for travelling. Insert a microSD card with 'Sound Oasis' sounds for unlimited sound library expansion. Has a Bluetooth wireless range of 30'/10m, & features built-in rechargeable lithium battery, allowing for 5 - 8 hours of cordless use on one charge.



Sound Oasis Bluetooth Speaker BST-100

| | | | | | |
|------------|---|--------|--------------------|--|--|
| Order code | 258-109 | Colour | Silver/Grey | | |
| Price | £44.80 Ex. Vat £53.76 Inc. Vat | Size | 63mmØx 66mm | | |

PERSONAL LISTENING & COMMUNICATION

LISTENING TO CONVERSATION

PERSONAL LISTENERS

If you have difficulty catching everything that is said in conversation you may benefit from using one of our Personal Listeners which amplify the sounds you want to hear and brings them closer to your ears.

HOW TO USE A PERSONAL LISTENER?

You can use a personal listener for conversation with friends or at work where you might find listening in a group or noisy situations stressful. Then, later in the day, you can use it to relax & listen to TV, radio, & music on your device or computer.

CHOOSE THE WAY YOU WANT TO LISTEN

Choose between headphones, headset or earphones, or use your hearing aid on the "T" Loop setting or Loop Program with a neckloop or ear hook. Your hearing aid may have the "T" (loop) facility. With this, you can enhance your listening by making more of your hearing aid. Sounds will be transmitted directly to your hearing aid, significantly reducing background noise.



SARABEC QUALITY

Sarabec® audio products provide clear, sharp sound reproduction leading to relaxed, stress free listening. These versatile products allow you to connect directly to various types of audio & visual equipment.




LISTENING OPTIONS

These listening options can be used with Crescendo, Sonumaxx & Earis. They can also be used independently with any piece of equipment that has a 3.5mm jack socket such as music players, laptops and tablets.

HEADSET

Designed to be worn under the chin. Lightweight and comfortable, these are the most popular headsets giving the loudest output. They are suitable for mild to moderate hearing loss. They have rubber ear tips that are very soft and comfortable to wear. These ear tips are washable and can be renewed.



HEADPHONES

Traditional style lightweight over the head headphones with an adjustable headband. Suitable for mild to moderate hearing loss. Can also be worn over hearing aids, and use either the microphone or "T" Loop setting or Loop Program.



EAR PHONES

Modern in-ear earphones designed to be lightweight, comfortable and discreet, they have the same mini loudspeakers as the headset. They are supplied with clear silicone ear moulds to help them stay in your ears and seal so that feedback can be reduced. Suitable for a mild to moderate hearing loss.



NECKLOOP

If your hearing aid has the "T" Loop setting or Loop Program, then you can make more of your hearing aid and its personalised settings by using a loop. The neckloop is a personal loop worn around the neck and transmits the sound directly to your hearing aid when switched to "T" or loop setting.



EAR HOOK

An ear hook or ear loop has a very small loop built-in to a hook similar to the shape of a behind the ear hearing aid. It is extremely thin and comfortable to wear. It is often found to be more discreet than a neckloop. Mono and Stereo versions are available. For those with two hearing aids, the stereo version would be used.





PERSONAL LISTENERS CRESCENDO 60 RANGE

WHAT IS THE CRESCENDO PERSONAL LISTENER?

The Crescendo 60 personal listeners will amplify conversation, or even the sound of television & audio, allowing you to listen at a volume that is comfortable for you. There is now no need to have the television or radio volume at a level that will disturb your family, friends or neighbours. You can take part in conversation with the help of the Crescendo which transmits the sounds clearly into your ears, effectively improving your

overall listening experience.

FOR USE WITH THE TV & FOR CONVERSATION



The Crescendo Personal Listener is a portable listener which is light and comfortable to use. For conversation, use the plug-in microphone to pick up the sound. For use with a TV or radio either a microphone or direct connection lead is used to pick up the sound. The Crescendo Personal Listener requires no technical knowledge to assemble, simply plug the accessories into the appropriate sockets and it is ready to use.

WHAT DO THE CRESCENDO SYSTEMS OFFER?

The range of systems listed allows you to choose a complete ready-to-use system, especially suited to your needs.

Some of the features of the Crescendo 60's include:

- A rotary volume control
- A rotary tone control (to boost high or low tones)
- An indicator "on" light with a low battery warning
- Mode selection microphone, direct input & loop pickup
- 2 AA alkaline batteries (supplied) for up to 60hrs use
- Belt clip & Neck cord
- A plug-in microphone for 1:1 communication
- Additional microphone (1.2m lead & tie-clip supplied) in case microphone needs to be placed closer to the sound source.
- Velcro pads for attaching the additional microphone to a loudspeaker grille
- A 3m extension for microphone or other listening option.

| Crescendo 60/1 system | | | |
|-----------------------|---|---|-------------------|
| Order code | 109-101 | Colour | White/Grey |
| Price | £79.99 Ex. Vat £95.99 Inc. Vat | Size | 95x62x25mm |
| | |   | |

Crescendo 60/1 system
Under-chin Stethoscope Headset, Plug-in Microphone, (1.2m lead), & 3m Microphone Extension Lead.



Telephone (voice): 01642 247789
Fax: 01642 230827 Textphone (NGTS): 18001 07497 242 294

www.sarabec.com
Email: mail@sarabec.co.uk

Crescendo 60/2 System

Supplied with Over-The-Head Headphones, Plug-in Microphone, Microphone (1.2m lead) & 3m Microphone Extension Lead.



Crescendo 60/2 system

Order code **109-102**
Price **£79.99 Ex. Vat**
£95.99 Inc. Vat

Colour **White/Grey**
Size **95x62x25mm**



Crescendo 60/5 System

Supplied with an NL6 Neckloop, Plug-in Microphone, Microphone (1.2m lead), 3m Microphone Extension Lead.



Crescendo 60/5 system

Order code **109-105**
Price **£79.99 Ex. Vat**
£95.99 Inc. Vat

Colour **White/Grey**
Size **95x62x25mm**



Crescendo 60/8 System

Supplied with an Under-chin Stethoscope Headset, NL6 Neckloop, Plug-in Microphone, Microphone (1.2m lead), 3m Microphone Extension Lead.



Crescendo 60/8 system

Order code **109-108**
Price **£95.99 Ex. Vat**
£115.19 Inc. Vat

Colour **White/Grey**
Size **95x62x25mm**



Crescendo 60 Personal Comms System

Includes Carry Case, Stethoscope Headset, 5 sets of spare ear tips, NL6 Neckloop, Plug-in Microphone, Microphone (1.2m lead), 3m Microphone Extension Lead and 2 AA batteries.



Crescendo Personal Comms System

Order code **109-107**
Price **£129.99 Ex. Vat**
£155.99 Inc. Vat

Colour **White/Grey**
Size **95x62x25mm**



röger

Roger is an intelligent wireless technology that transmits speech directly to hearing aids and cochlear implants, helping to overcome distance and noise.

Functioning like a third ear, Roger allows people with hearing loss to learn, work and socialise more easily - a big step towards more inclusive societies and greater personal well-being. Roger technology works as a system that enables a speaker's voice to be transmitted from a Roger wireless microphone to a listener's hearing aid or cochlear implant via Roger receivers.

Up to **62%** more understanding of speech in noise and over distance when wearing hearing aids in combination with Roger microphones, compared to people without hearing loss.



Phonak Roger On Transmitter (Champagne)

Phonak's latest state-of-the-art wireless microphone **Roger On** incorporates MultiBeam 2.0 Technology, providing excellent speech-understanding in noisy situations. Roger On is a versatile microphone transmitter, ideal for situations where background noise is present such as pubs and restaurants. It is used in conjunction with suitable Roger receivers, such as the Roger Neckloop receiver or Roger X ear-level receivers (purchased separately).



Phonak Roger On Transmitter (Champagne)

| | | | |
|------------|---|--------|--------------------|
| Order code | 240-004 | Colour | Champagne |
| Price | £669.00 Ex. Vat £802.80 Inc. Vat | Size | 100x24x14mm |



Phonak Roger Neckloop (03)

Phonak Roger Neckloop is a universal inductive neckloop receiver that allows any person with hearing loss to have industry-leading speech-in-noise benefits.

Phonak Roger Neckloop (03) can be used together with all Roger microphones, and it is compatible with hearing aids, BAHA's, and cochlear implants with a "T" Loop Setting or Program, regardless of the hearing device brand. It is worn around the neck by adults, students, and children above 8 years old to give access to the complete benefits of Roger technology in noisy situations and over distance, bringing conversations closer.

Designed for everyday use – with its lightweight and detachable neckloop in two lengths, it is comfortable to use. Once connected to a Roger microphone it stays connected. Just switch it on and it lasts for 10 hours. Roger Neckloop uses a USB C interface for charging and digital output, has an easy-to-read display, and a 3.5mm headphone audio output.



Phonak Roger Neckloop (03)

Roger Neckloop replaces Roger Mylink.

| | | | |
|------------|---|--------|-------------------------|
| Order code | 240-048 | Colour | Silver/Champagne |
| Price | £134.50 Ex. Vat £161.40 Inc. Vat | Size | 66x29x12mm |
| | | | |

Roger Select (Champagne)

The Roger Select is a versatile, wireless transmitting microphone from Phonak that uses their innovative MultiBeam Technology. The Select has multiple microphones pointing in six directions covering speech over 360°.

When it is placed in the centre of a table, it automatically selects the person who is talking and seamlessly switches from one speaker to another, keeping you part of the conversation. When there are individual conversations around the table you can tap the segment in the direction of the person you are speaking with to focus on their speech.

Roger Select can be used in many different situations such as:

- In restaurants, cafés, etc.
- Small group conversations with up to 6 participants
- 1:1 conversations with a friend, partner or caregiver
- During lectures or in class
- Outdoor and indoor activities (gym, guided tour)
- Watching TV
- Landline phonecalls (with adaptor - not included)
- Audio/voice from the computer, radio



Roger Select (Champagne)

| | | | |
|------------|---|--------|------------------|
| Order code | 240-009 | Colour | Champagne |
| Price | £624.00 Ex. Vat £748.80 Inc. Vat | Size | 55x12mm |
| | | | |

Phonak TV Connector

The Phonak TV Connector is a wireless hearing aid accessory that streams TV and music directly to your hearing aids, transforming your Bluetooth enabled hearing aids into wireless stereo headsets. A simple plug-and-play interface works with TVs and other audio sources to support direct hearing aid connectivity to Lumity, Paradise, and Marvel hearing aids.

Phonak TV Connector

| | | | |
|------------|---|--------|-------------------------|
| Order code | 240-127 | Colour | Black |
| Price | £159.00 Ex. Vat £190.80 Inc. Vat | Size | 50mmx 50mm x15mm |
| | | | |



PERSONAL AUDIO ACCESSORIES

Explore our range of Listening accessories for use with the systems we supply, and for those that have now been discontinued.

AMPLIFY SOUND WITH LISTENING ACCESSORIES

You may not need an amplifying system. You might find that you can improve your entertainment enjoyment simply by utilising one of these accessories, and plugging them directly into your own system. Typically, you can use any of these accessories where there is a suitable headphone output. Although connecting the accessories to your own equipment as opposed to one of our full systems may not provide the same level of power, it may still become loud enough for you to use more comfortably.

FOR USE WITH OR WITHOUT HEARING AIDS

Using any of these products with or without a hearing aid will help bring the sound closer to you, while reducing background noise. For hearing aid users, the sound is transmitted inductively, directly to the hearing aid and thereby cuts out background noise.

USING YOUR HEARING AIDS LOOP PROGRAM

If you want to use your hearing aids "T" Loop setting or Loop Program, then we would recommend the Ear Hooks as these give excellent reproduction. If you have two hearing aids and use the stereo Ear Hook, then please note that if your MP3 player has stereo output you will hear in Stereo, unlike a Neckloop which only transmits in Mono. For these items, headphone extension leads can be used to extend the distance between you and the audio equipment.



HEADPHONES

LH-056 Headphones

The LH-056 headphones have integrated volume control in the connecting lead, which allows the separate adjustment of the left and right side. This lightweight headphone is specially designed for speech intelligibility. They are very comfortable, and will give a relaxing sound experience. The connecting lead is 5.5m long and comes with a 3.5mm jack plug and a 6.3mm adapter.

| LH-056 Headphones | | | |  |
|-------------------|------------------------|--------|--------------|---|
| Order code | 230-296 | Colour | Black |   |
| Price | £47.39 Inc. Vat | | | |

HS-101 Stereo Headphones

A pair of hearing aid-compatible traditional style, lightweight, over-the-head headphones with an adjustable headband. Suitable for mild to moderate hearing loss. These headphones contain a telecoil so the user can change their hearing aid to the "T" Loop setting or Loop Program, and the sound will be transmitted directly to the hearing aids.

- 3.5mm Stereo Jack plug
- 1.2m Lead length
- 32 Ohm Impedance



| HS-101 Stereo Headphones | | | |  |
|--------------------------|------------------------|--------|--------------|---|
| Order code | 303-061 | Colour | Black |   |
| Price | £16.79 Inc. Vat | | | |

HEADSETS

Stethoscope Headset HS7

Designed to be worn under the chin. Lightweight and comfortable, these are the most popular headsets giving the loudest output. They are suitable for mild to moderate hearing loss. They have rubber ear tips that are very soft and comfortable to wear. These ear tips are washable and can be renewed.

- 1.0m Connecting Lead
- 3.5mm Stereo Jack Plug
- Impedance:16 ohm

| Stethoscope Headset HS7 | | | |  |
|-------------------------|------------------------|--------|--------------|---|
| Order code | 307-107 | Colour | Black |  |
| Price | £14.99 Inc. Vat | | | |



NECKLOOPS

Neckloop NL6 - 75cm Stereo Lead

The Neckloop is designed for use with our Earis, Sonumaxx, Crescendo 60, 50 and 20 Plus systems. The Neckloop is for use with hearing aids that have the "T" Loop setting or Loop Program, and can be used with most TVs, radios or other equipment that have a 3.5mm stereo headphone socket.

Not suitable for use with the Crescendo 10 or 20.

| Neckloop NL6 | | Colour Black | | |
|--------------|------------------------|---------------------|--|--|
| Order code | 201-126 | | | |
| Price | £22.49 Ex. Vat | | | |
| | £26.99 Inc. Vat | | | |

Geemarc iLoop+ Neckloop with Microphone

The Geemarc CL iLoop+ is a powered neckloop that connects to any smartphone, iPhone, MP3, Kindle, etc. that has a 3.5mm headset socket. Using 2 x AAA rechargeable batteries (included), it can provide up to 30dB of amplification. Switch between music and phonecalls at the touch of a button.

| Geemarc iLoop+ Neckloop | | Colour Black | | |
|-------------------------|------------------------|---------------------|--|--|
| Order code | 252-187 | | | |
| Price | £54.16 Ex. Vat | | | |
| | £64.99 Inc. Vat | | | |

EAR HOOKS

Music Link - Stereo/Mono

The Music Link, otherwise known as an ear hook or ear loop, is a miniature telecoil inside a very thin and comfortable hook that sits behind the ear, but can equally be used with an in-the-ear hearing aid with the "T" Loop setting or Loop Program. Use the mono for one hearing aid and the stereo for two hearing aids.

| Music Link - Mono Ear Hook | | Colour Black/Silver | | |
|----------------------------|------------------------|----------------------------|--|--|
| Order code | 315-001 | | | |
| Price | £19.99 Ex. Vat | | | |
| | £23.99 Inc. Vat | | | |

| Music Link - Stereo Ear Hook | | Colour Black/Silver | | |
|------------------------------|------------------------|----------------------------|--|--|
| Order code | 315-002 | | | |
| Price | £29.99 Ex. Vat | | | |
| | £35.99 Inc. Vat | | | |

MICROPHONES

Standard Microphone 1.2m lead

Standard microphone for many of our products including Crescendo and LA loop amplifiers.

- Electret Microphone
- 3.5mm jack plug
- 1.2m lead length

| Standard Microphone 1.2m lead | | Colour Black | | |
|-------------------------------|------------------------|---------------------|--|--|
| Order code | 330-005 | | | |
| Price | £17.99 Inc. Vat | | | |
| | | | | |

Plug In / Communicator Microphone

Communicator Microphone, which gives the same performance as our standard 1.2m microphone, but doesn't have a trailing lead. Sometimes described as a "Stubby or Plug-In" microphone.

- 3.5mm Audio Jack Plug
- Electret microphone

| Plug In / Communicator Microphone | | Colour Black | | |
|-----------------------------------|------------------------|---------------------|--|--|
| Order code | 330-009 | | | |
| Price | £13.79 Inc. Vat | | | |
| | | | | |

HOME ALERTING

Alerting - Telling you what's happening

WHAT ARE ALERTING DEVICES?

Alerting devices make you aware of everyday occurrences at home or at work, e.g. the ringing of the telephone or doorbell, alarm clocks, the sound of the smoke or fire alarm, or a baby crying. You'll either "feel" these occurrences via a vibrating pad or "see" that these are happening via a flashing light. This Alerting section offers a rich variety of choice.

THE COMPLETE HOME ALERTING SYSTEM

Full home alerting systems allow you to completely tailor a system to your needs - meaning a single receiver device will warn you of all the occurrences mentioned above. Add more receivers to these systems for added flexibility.

SMOKE ALARM SYSTEMS

Smoke alarm systems are invaluable to hard of hearing people, particularly if they remove their hearing aids at night. You could benefit from the independence and security of such a system with a vibrating pad and flashing strobe.

DOORBELL SYSTEMS

A range of totally wire-free doorbells, either battery or mains powered, means that you

won't miss that important parcel or caller at your door. You can use the portable chimes in whichever room you choose, and the battery powered chimes can be taken into the garden or shed.

Please remember that building structures e.g. metal, uPVC door frames and thick walls may affect transmission. **Where a range for a product is quoted this is the maximum without obstructions.*

TELEPHONE SYSTEMS

Do you miss out on telephone calls because you don't hear the phone ringing? No matter if you're at work in a noisy workshop, watching TV in your lounge with the doors closed and the phone is ringing in the hallway, or you're listening to music through your hearing aid or headphones, these telephone alerting products can help keep you alerted.

ALARM CLOCKS

Getting up in the morning can be a chore for anyone, but if you have difficulty hearing the alarm, your sleep will be less than peaceful if you are worried about sleeping in. A little assistance from either a vibrating pad and/or a flashing light could make that morning feel less stressful.

Telephone (voice): 01642 247789
Fax: 01642 230827 Textphone (NGTS): 18001 07497 242 294

www.sarabec.com
Email: mail@sarabec.co.uk



THE SIGNOLUX ALERTING RANGE

Get alerted to events around the home with our complete Signolux systems. Explore our range of transmitters, receivers, accessories & batteries.

The Signolux alerting systems now include additional features found on more expensive systems to make it the affordable solution for deaf and hard of hearing people alike. The simple, efficient, and cost effective system is ideal for those who are wanting peace of mind.

Get alerted to sounds around the home. The receivers will alert you via a combination of loud sound, visual indications and or vibrations. This sophisticated system uses just two universal transmitters to keep installation as simple as possible. The options include an audio transmitter that uses microphones to detect sound, or a direct transmitter that can be connected to existing devices around the house. *Dedicated transmitters are available for person-to-person, bell push, and smoke alarm.*



FEATURES AND BENEFITS:

- All receivers come with either a vibrating pad, strobe for bright white flashing lights and coloured indicators, or loud sound alerts.
- Receivers have either 6 or 8 coloured icons to visually identify the source of the sound
- Up to 8 transmitters can be coded to each icon
- Choice of 15 loud sounds up to 90dB output
- Wirefree; uses 868 MHz radio frequency to send and receive system signals
- Secure coding system to avoid interference from neighbouring systems
- Operating range of up to 200m
- Wireless receivers can be carried around for use in and around the house
- Portable pocket or clip-on vibrating receiver
- Plug-in mains powered strobe receiver
- Easy wireless installation between transmitters (sound detectors) and receivers
- Additional receivers and transmitters can be



- added to your system as required
- Battery operated with mains adaptor option
- Batteries included with the kits where necessary



SIGNOLUX SYSTEMS

Signolux Tower Doorbell System

Introducing the new standard system in the Signolux range, designed for people with hearing problems to inform about domestic events. The 2-part system, made up of a Signolux bellpush and Signolux Tower receiver, comes supplied with rechargeable batteries and a mains power supply as standard.

- 16 x chime sounds
- Extra loud: 95dB+ on some sounds
- Mains powered + rechargeable
- 3.5mm socket for vibrating pad & other alerts
- 8 x easy-to-see visual indicators
- Bright flashing light with 360° visibility
- Compatible with older Signolux transmitters



Signolux Plug-in Doorbell System

A plug-in flashing doorbell system which also interacts with the standard Signolux system transmitters.

- Indicators - flashing lights & coloured LEDs
- 3.5mm socket for vibrating pad & other alerts
- Doorbell symbol displayed on the front
- Displays telephone, baby cry, pager, fire & alarm alerts when attached to the correct transmitters
- Supports additional signal transmitter modules



Please note, this is a flashing/indication ONLY system, it does not make any sound unless a sound module is plugged into it.

Signolux Tower Doorbell System

Includes: Flash & Chime Tower Receiver, Bellpush, Rechargeable Batteries & Charger.

| | | | |
|------------|---|--------|-----------------|
| Order code | 230-750 | Colour | White |
| Price | £72.49 Ex. Vat £86.99 Inc. Vat | Size | 143x80mm |



Signolux Plug-In Doorbell System

Includes: Plug-in Flash Receiver, Bellpush, & Batteries.

| | | | |
|------------|--|--------|---|
| Order code | 230-714 | Colour | White |
| Price | £94.33 Ex. Vat £113.20 Inc. Vat | Size | Bellpush: 80x40x16mm Receiver: 120x65x45mm |



Signolux Tower System & Vibrating Pad

The Signolux Tower Doorbell System + Vibrating Pad was designed for people with hearing problems and informs about domestic events with a bright 360° flashing light, LED display of the event symbols, & powerful tones. The vibrating pad slips under a pillow or mattress, and alerts via strong vibration patterns.



Signolux Tower System & Phone Transmitter

Get alerted to domestic events via bright 360° flashing lights, LED display of the event symbols, and powerful tones.



Signolux Tower System & Vibrating Pad

Includes: Signolux Tower Receiver, Bellpush, Vibrating Pad, Rechargeable Batteries & Charger.

| | | | |
|------------|--|----------------|-----------------|
| Order code | 230-751 | Colour | White |
| Price | £90.49 Ex. Vat £108.59 Inc. Vat | Receiver Size: | 143x80mm |
| | | Vibrating Pad: | 90x28mm |



Signolux Tower System & Phone Transmitter

Includes: Signolux Tower System (230-750), Phone Transmitter, Rechargeable Batteries, and Charger.

| | | | |
|-------|---|--------------|--------------------|
| Code | 230-752 | Colour | White |
| Price | £108.49 Ex. Vat £130.19 Inc. Vat | Receiver: | 143x80mm |
| | | Transmitter: | 122x75x38mm |



COMING SOON!

connect

INTRODUCING
THE BRAND NEW
SIGNOLUX CONNECT!
COMING SOON.

With the **Signolux Connect T**, you'll always be aware of events happening in the home. Simply connect your smartphone to the Connect T via Wi-Fi or Bluetooth, and never miss an alert again.



SIGNOLUX TRANSMITTERS

Signolux Audio Universal Transmitter

Used to pick up audio signals from different sources such as a baby crying, an existing doorbell, an intercom system or the telephone ringing.



Signolux Direct/Phone Universal Transmitter

Used to alert the user to the telephone ringing (by plugging into a BT socket), can be used to wire into an existing hard-wired doorbell system, alarm panel or water level indication.



Signolux Audio Universal Transmitter

| | | | |
|------------|---|--------|--------------------|
| Order code | 230-720 | Colour | White |
| Price | £34.52 Ex. Vat £41.42 Inc. Vat | Size | 122x75x38mm |



Signolux Direct/Phone Universal Transmitter

| | | | |
|------------|---|--------|--------------------|
| Order code | 230-721 | Colour | White |
| Price | £34.52 Ex. Vat £41.42 Inc. Vat | Size | 122x75x38mm |



Signolux Sound Monitor Transmitter

The Signolux Sound Monitor is capable of learning 3 different sounds such as tones from timers, ringtones, alarm, sounds from CO2-detectors & smoke alarms etc. Once learned, the system then transmits these sounds as alerts via a Signolux receiver.



Signolux Baby Monitor Audio Transmitter

Can be used to pick up audio signals from a baby crying, and has a sensitivity control, so you can set it to the strength of your baby's cry. This transmitter sends a radio signal to a Signolux receiver which will alert you through audio and visual signals.



Signolux Sound Monitor Transmitter

Includes: Rechargeable Batteries, and Charger

| | | | |
|------------|---|--------|--------------------|
| Order code | 230-718 | Colour | White |
| Price | £119.90 Ex. Vat £143.88 Inc. Vat | Size | 122x75x38mm |



Signolux Baby Monitor Audio Transmitter

Includes: Plug-in microphone and Charger

| | | | |
|------------|---|--------|--------------------|
| Order code | 230-719 | Colour | White |
| Price | £43.90 Ex. Vat £52.68 Inc. Vat | Size | 122x75x38mm |



SIGNOLUX BELLPUSH

Signolux Bell Push / Person Call Button

A vertical door bell push for the Signolux system which can be used at either the front or back door, or works well as a personal call button. The push features a red confidence light and can be wall-mounted using the supplied double-sided adhesive fixing strip.



Signolux Bellpush Landscape/ Person Call Button

The Signolux wireless landscape Bell Push can be used as an additional Doorbell Push for a Signolux Visual Alert System or as a replacement.

- Weather resistant to IP55
- LED Confidence Light
- Name Plate



Signolux Bell Push / Person Call Button

Bell push with CR2032 battery

| | | | |
|------------|---|--------|-------------------|
| Order code | 230-710 | Colour | White/Grey |
| Price | £16.15 Ex. Vat £19.38 Inc. Vat | Size | 80x40x16mm |



Signolux Bellpush Landscape/ Person Call Button

Landscape bell push with CR2032 battery and name card

| | | | |
|------------|---|--------|-------------------|
| Order code | 230-711 | Colour | Grey/White |
| Price | £16.15 Ex. Vat £19.38 Inc. Vat | Size | 80x40x16mm |



SIGNOLUX DETECTORS

Signolux Guardion Smoke Detector

The Signolux Guardion Smoke Detector works seamlessly with the Signolux Visual Indication System. It operates as a standard smoke alarm, but also sends a wireless signal to a Signolux receiver so that a user may be alerted by loud sounds, visual indications and vibrations.

- Optical smoke alarm transmitter
- Fitted (non-removable) lithium battery



Signolux Guardion Smoke Detector

| | | | |
|------------|--|--------|-------------------|
| Order code | 230-724 | Colour | White |
| Price | £89.99 Ex. Vat £107.99 Inc. Vat | Size | 95mmØx50mm |



Signolux Guardion CO/Heat Detector

This combined CO and heat detector detects and reports elevated carbon monoxide levels in the room air at an early stage and also reacts when the room temperature exceeds a limit of 60 °C. Particularly great for rooms with fireplaces or garages or in areas in which there are fuel-operated devices and incomplete combustion processes which cannot be ruled out, or where there is a risk that harmful smoke gases like carbon monoxide can spread. Connects to any Signolux receiver in the Signolux range. Ideal for use in areas unsuitable for regular smoke detectors as they feature optoelectronic sensor technology (smoke particle detection), which can increasingly lead to false alarms. Up to 10-year battery life (fixed battery, not changeable).

- Self-testing function & warning given at low battery level
- Loud alarm tone



Signolux Guardion CO/Heat Detector

| | | | |
|------------|--|--------|-------------------|
| Order code | 230-728 | Colour | White |
| Price | £89.99 Ex. Vat £107.99 Inc. Vat | Size | 95mmØx50mm |



Signolux Alarmo 2 Detector

When an independent smoke alarm sound is detected, Alarmo transmits a radio signal to all Signolux receivers within range, which then alerts you to the event.

There is no need to replace your smoke alarm because Alarmo can detect the alarm sounds from existing smoke alarms without the need for any modification or wiring.

It incorporates a highly intelligent processor which listens out for a typical domestic smoke alarm sound, and is therefore extremely good at not creating false alerts.



Signolux Alarmo 2 Detector

| | | | |
|------------|---|--------|--------------------|
| Order code | 230-727 | Colour | White |
| Price | £119.90 Ex. Vat £143.88 Inc. Vat | Size | 112x95x38mm |



SIGNOLUX RECEIVERS



Signalux Tower Receiver

Introducing the new standard receiver in the Signalux range. It comes supplied with rechargeable batteries and a mains power supply as standard. This product was designed for people with hearing problems and informs about domestic events with a bright 360 degree flashing light, LED display of the event symbols and powerful tones, the volume of which can be adjusted.

- 16 x chime sounds
- Extra loud: 95dB+ on some sounds.
- Mains powered + rechargeable
- 8 x easy-to-see visual indicators
- 3.5mm socket for vibrating pad & other alerts
- Bright flashing light with 360° visibility
- Compatible with older Signalux transmitters.

Signalux Tower Receiver

| | | | |
|------------|---|--------|-----------------|
| Order code | 230-713 | Colour | White |
| Price | £56.49 Ex. Vat £67.79 Inc. Vat | Size | 143x80mm |



Signalux DS2 Alarm Clock

The DS2 drives all of your alerts to a single point; your bedside table. The device receives all types of signals from the system and alerts to whatever is going on. The DS2 provides connectors so you can connect to any of the Signalux receiver alerts. It offers a single piece of equipment for all of your alerting needs; telephone, door, smoke and alarm clock. It provides night time coverage and will alert you to any alerts in the home. It utilises intelligent digital technology, offering a wide range of waking methods. A particularly powerful sound generator with an exceptionally bright row of built-in flashing LED's.

- Large, easily discernible illuminated digits for time and alarm time display
- Unmissable acoustic alarm signal
- The alarm tone pitch and volume are adjustable
- Visual alarm with a row of flashing LEDs
- Snooze-function
- Display brightness adjusts automatically to its

environment

- 12 or 24 hour time display
- Back-up against power failure (4 x AAA batteries not supplied)
- 3.5mm socket for vibrating pad & other alerts



Signalux DS2 Alarm Clock

| | | | |
|------------|---|--------|---------------------|
| Order code | 230-708 | Colour | White |
| Price | £159.90 Ex. Vat £191.88 Inc. Vat | Size | 123x210x63mm |



Signalux Plug-in Flash Receiver

Plug-in mains receiver for use in kitchens or for any mains socket that is placed at eye level.

- Bright white flashing light
- 6 coloured visual notification icons
- 3.5mm socket for additional alerts

Signalux Plug-in Flash Receiver

| | | | |
|------------|---|--------|--------------------|
| Order code | 230-705 | Colour | White |
| Price | £80.55 Ex. Vat £96.66 Inc. Vat | Size | 120x65x45mm |





Signolux Portable Vibrating Receiver

The Signolux Portable Receiver is a portable vibrating pager unit which will alert you through both lights and vibrations. It can be carried anywhere in your home and can also be used at night time when used in conjunction with an external vibrating pad (sold separately).

- Belt clip and neck strap
- Receiver vibrates and displays flashing icons upon signal being received
- 6 coloured visual notification icons
- 3.5mm socket for vibrating pad & other alerts
- Optional charger with rechargeable batteries



Signolux Portable Vibrating Receiver

| | | | | |
|------------|------------------------|--------|-------------------|--|
| Order code | 230-706 | Colour | White | |
| Price | £80.55 Ex. Vat | Size | 90x58x35mm | |
| | £96.66 Inc. Vat | | | |

Signolux Portable Vibrating Receiver Pack

The Signolux Portable Receiver is a portable vibrating pager unit which will alert you through both lights and vibrations. It can be carried anywhere in your home and can also be used at night time when used in conjunction with the charger and external vibrating pad.

The charger base station for the Signolux Portable Vibrating Receiver is supplied with rechargeable batteries. The lead from the vibrating pad is plugged into the side of the pager receiver and the pad is placed under your pillow, ready to vibrate when the alarm goes off. The vibrating pad is powered by 3 x AA alkaline batteries (supplied). This device also has a low battery indicator. The pack includes:

- Vibrating Pager Receiver
- Pager Charger
- Vibrating Pad
- Power supply



Signolux Portable Vibrating Receiver Pack

| | | | | |
|------------|-------------------------|--------|--|--|
| Order code | 230-715 | Colour | White | |
| Price | £135.90 Ex. Vat | Size | Portable Receiver: 90x58x30mm | |
| | £163.08 Inc. Vat | | Charger: 75x80x65mm Vibrating Pad: 82x82x30mm | |

Signolux Portable Vibrating Receiver Charger

Charger base station for the Signolux Portable Vibrating Receiver. Comes supplied with:

- Rechargeable batteries
- 2 x 3.5mm sockets for vibrating pad & other alerts
- Mains power supply

Signolux Portable Vibrating Receiver Charger

| | | | | |
|------------|------------------------|--------|-------------------|--|
| Order code | 230-707 | Colour | White | |
| Price | £39.90 Ex. Vat | Size | 75x80x65mm | |
| | £47.88 Inc. Vat | | | |



Home Alerting

SIGNOLUX ACCESSORIES

Vibrating Pillow c/w Batteries

This Vibrating Pad can be connected to any of the Signolux receivers. This equipment uses its internal 3 x AA batteries so does not draw power from the connected receiver, which is especially useful if using the system in a portable environment. 1.5m lead.



Signolux Vibrating Pad (un-powered)

An un-powered vibrating pad for use with any of the Signolux receivers (except 230-705). 1.5m lead.



Vibrating Pillow c/w batteries

| | | | |
|------------|------------------------|---|-------------------|
| Order code | 230-315 | Colour | Black |
| Price | £18.49 Ex. Vat | Size | 82x82x30mm |
| | £22.19 Inc. Vat |   | |

Signolux Vibrating Pad (un-powered)

| | | | |
|------------|------------------------|---|------------------|
| Order code | 230-735 | Colour | Black |
| Price | £15.15 Ex. Vat | Size | 90Øx 28mm |
| | £18.18 Inc. Vat |   | |

Flash Module MF-1

Emits signals via powerful flashes of light; they make the whole room brighter and can be detected even with your eyes closed. Can be connected to any receiver from the Signolux range.



Sound Module MA-1

Add to any existing Signolux receiver to hear the alerts as well as see and feel them. Generates a powerful alarm tone - continuously variable tone and volume settings.



Flash Module MF-1

| | | | |
|------------|------------------------|---|---------------------|
| Order code | 230-317 | Colour | Silver |
| Price | £56.17 Ex. Vat | Size | 65x93x 134mm |
| | £67.40 Inc. Vat |   | |

Sound Module MA-1

| | | | |
|------------|------------------------|---|---------------------|
| Order code | 230-318 | Colour | Silver |
| Price | £49.81 Ex. Vat | Size | 134x93x 65mm |
| | £59.77 Inc. Vat |   | |

Signolux Alarm Pressure Mat

An Alarm Pressure Mat that can be used with the Signolux Direct Universal Transmitter. It is a durable pressure pad sensor that can be used for various applications (floor pressure mat, bed leaving, chair leaving etc.) 3m Lead.



Signolux Water Level Sensor

Water Level Sensor for use with Signolux Direct Transmitter. Suitable for almost any liquid. Attach the sensor to the edge of e.g. a wall using the double-sided sticky tape (included), so that the contacts are touching the floor. The connected transmitter emits a signal if liquid touches both of the contacts. 2m lead.



Signolux Alarm Pressure Mat

| | | | |
|------------|------------------------|---|------------------|
| Order code | 230-731 | Colour | White |
| Price | £40.27 Ex. Vat | Size | 750x450mm |
| | £48.32 Inc. Vat |   | |

Signolux Water Level Sensor

| | | | |
|------------|------------------------|---|--------------|
| Order code | 230-732 | Colour | Black |
| Price | £38.90 Ex. Vat |   | |
| | £46.68 Inc. Vat | | |

Contact Switch (Normally Open)

To use with Signolux Direct Transmitter for detecting door or window opening.





Signolux Power Supply (7.5v 1.0A)


Mains Power supply for Signolux system transmitters and receivers. Replaces part No. 230-722



Contact Switch (Normally Open)

| | | | |
|------------|------------------------|---|--------------|
| Order code | 230-695 | Colour | White |
| Price | £14.40 Inc. Vat |   | |

Signolux Power Supply (7.5v 1.0A)

| | | | |
|------------|------------------------|---|--------------|
| Order code | 230-792 | Colour | Black |
| Price | £10.16 Inc. Vat |   | |



SIGNOLUX TOWER BABY MONITOR SYSTEM



£119.90
Ex. Vat

GET ALERTED TO YOUR BABY'S CRY FROM ANYWHERE IN THE HOME.

Comprised of: Signolux Tower Receiver, Vibrating Pad, Signolux Audio Transmitter, Plug-In Microphone, & Mains Power Supply.

- Portable for use throughout the home
- Rechargeable & mains powered/charged
- Operating range - 200m on the 868Mhz frequency.
- 8 indicators meaning that you can be alerted to different events within the home
- 16-tone sequences/melodies for acoustic alerting



| Signolux Tower Baby Monitor System | |
|---|--|
| Includes: Flash & Chime Tower Receiver, Bellpush, Rechargeable Batteries & Charger. | |
| Order code | 230-753 |
| Price | £119.90 Ex. Vat £143.88 Inc. Vat |
| Colour | Black/White |
| Size | Receiver: 143mm x 80mm Transmitter: 122x75x38mm |
| | |



Wi-Safe® 2

HOME SAFETY SOLUTIONS FOR THE HEARING IMPAIRED

WHAT IS THE WI-SAFE2 RANGE?

The new generation of FireAngel Wi-Safe2 interlinked safety products are designed to provide the ultimate in-early warning to escape the dangers of fire and carbon monoxide. They are the ultimate cost-effective solution when it comes to protecting vulnerable groups or homeowners with specific needs from the dangers of fire and carbon monoxide (CO) poisoning. Wi-Safe2 ensures homeowners feel safe from danger and can live independently within their home to improve their quality of life.

HOW DOES THE THE WI-SAFE2 RANGE WORK?

Wi-Safe2 range incorporates heat, smoke, and carbon monoxide alarm detectors, which send wireless signals, not only amongst themselves, but to any combination of Flashing Strobe with Vibrating Pad and / or Loud Low Frequency Sounder. These will alert the homeowner to a potential fire or carbon monoxide incident.

WHO CAN BENEFIT FROM THIS PRODUCT?

The Wi-Safe2 products can be tailored to suit

individual needs thus ideal for everyone, but particularly for people with hearing impairments, and the disabled who may have reduced capacity to be alerted to danger within their home and escape in time.

WIRELESS AND BATTERY-POWERED

Due to their wireless operation, FireAngel's Wi-Safe2 products do not require costly or disruptive wiring in the home. Alarms are battery-operated, and so they are quick and easy to install.

FEATURES:

- Wireless products are quick and easy to install
- Unique coding system to avoid any false alarms from neighbouring properties
- The range can be tailored to suit needs
- Provides early warning to fire/carbon monoxide
- Improves the quality of life of homeowners
- Intelligent technology - creates a network of all the installed units in the home, so they all communicate with each other at all times.
- When ANY alarm detects an emergency or sounds - ALL alarms and receivers will sound and indicate an emergency.
- Intelligent 'locate' feature - on activation, press the test/reset button on any alarm/receiver to silence all but the initiating alarm.

Should a failure of your network occur, then a fault condition will be shown on either the strobe or loud sounder receivers.

RECOMMENDED PRODUCT POSITIONING



SMOKE ALARM

Landing / Hallway, Bedroom, Living Room
Minimum of one alarm should be installed on every floor of your home.




HEAT ALARM

Kitchen, Garage, Loft



CO ALARM

Kitchen, Living Room, Bedroom, Boiler Room
One alarm should be installed in every room containing a fuel-burning appliance.



Wi-Safe2 Strobe & Vibrating Pad Receiver

The combined strobe and pad receiver will respond to all Wi-Safe2 alarms in your network. The strobe will flash and the pad vibrate. The red LED illuminates immediately on activation of an alarm and shows whether the danger is either from fire or CO. The pad can be placed under a pillow or cushion and the strobe laid flat or wall-mounted. The strobe main unit button can be used to test all the alarms in your network - Great for anyone struggling to reach the ceiling. Supplied with mains power supply and a backup battery.



Wi-Safe2 Strobe & Vibrating Pad Receiver W2SVP630

| | | | |
|------------|-----------------|--------------|------------|
| Order code | 256-110 | Colour | White |
| Price | £91.26 Inc. Vat | Strobe size: | 102x150x40 |
| | | Pad: | 98 x 34 |



Wi-Safe2 Wireless Thermoptek Smoke Alarm WST-630

Battery powered, sealed-for-life, wireless interlink smoke alarm with Thermoptek technology. The WST-630 has diagnostics functionality which monitors alarm activations and facilitates fault-finding. Combines optical sensing with thermal enhancement to react to both fast flaming and slow smouldering fires, and so removes the need for Ionisation alarms that can have disposal issues.



Wi-Safe2 Thermoptek Smoke Alarm WST-630

| | | | |
|------------|-----------------|--------|----------|
| Order code | 256-174 | Colour | White |
| Price | £55.38 Inc. Vat | Size | 132 x 41 |



Wi-Safe2 Wireless Thermistek Heat Alarm WHT-630

Battery powered, sealed-for-life, wireless interlink heat alarm with Thermistek and built-in diagnostics technology. The WHT-630 is designed for areas where dust and fumes may trigger nuisance alarms with conventional smoke alarms - ideal for kitchens, garages, and lofts. Offers additional protection when used in conjunction with WST-630 smoke alarms.



Wi-Safe2 Thermistek Heat Alarm WHT-630

| | | | |
|------------|-----------------|--------|----------|
| Order code | 256-176 | Colour | White |
| Price | £55.31 Inc. Vat | Size | 132 x 60 |



Wi-Safe2 Wireless Carbon Monoxide Alarm W2CO10XT

Battery powered, sealed-for-life, wireless interlink carbon monoxide (CO) alarm with built-in diagnostics technology gives protection from CO poisoning in the home. The W2-CO-10XT contains advanced sensors to accurately detect levels of CO to provide early warning and better chance of escape.



Wi-Safe2 Carbon Monoxide Alarm W2CO10XT

| | | | |
|------------|-----------------|--------|---------------|
| Order code | 256-177 | Colour | White |
| Price | £62.10 Inc. Vat | Size | 114 x 95 x 40 |



FireAngel Digital CO Detector with Digital Display

The FireAngel FA3322 is a feature-rich carbon monoxide alarm with an advanced digital LCD display, powered by a tamper-proof 10 year lithium battery. An intelligent sensing algorithm is utilised, which increases the sensing rate when a low level of CO is detected, providing a faster response. The alarm's smart multi-function LCD display provides visual instructions for the user. The display illuminates and changes colours depending on the event type (e.g. standby, alarm etc). The FA3322 also includes environmental monitoring features.



FireAngel Digital CO Detector with Digital Display

| | | | |
|------------|-----------------|--------|---------------|
| Order code | 256-190 | Colour | White |
| Price | £39.24 Inc. Vat | Size | 129 x 76 x 31 |





DOORBELLS & ACCESSORIES

WIREFREE DOORBELLS

Never miss a visitor with our wide range of Chimeflash, Honeywell, and Signolux doorbell systems. Explore our totally wirefree and wireless doorbells, either battery or mains powered, and ensure you never miss that important parcel or caller at your door again. You can use the portable chimes in whichever room you want, and you can move them from room to room as they are portable and built for your convenience.

HASSLE-FREE

There is no need to worry if your neighbour has an identical system, as our doorbells come with a simple coding system to ensure that it is only your system that rings.

TRANSMISSION

Please remember that building structures e.g. metal, uPVC door frames and thick walls may affect transmission.

THE RANGE

Where a range for a product is quoted - this is the maximum without obstructions.

CHIMEFLASH DOORBELL & ACCESSORIES

SARABEC CHIMEFLASH DOORBELL

EXTRA LOUD AND ULTRA-BRIGHT ALERTING DOORBELL, PERFECT FOR d/DEAF, HARD OF HEARING, AND VISUALLY IMPAIRED PEOPLE.

This wireless, extremely loud doorbell system offers ultra-bright flashing LED lights, volume control, & indicator symbols. Perfect for d/Deaf, hard of hearing, & visually impaired individuals, or anyone wanting a loud doorbell with excellent range & reliability!

- Extra loud: Up to 100dB
- Up to three transmitters from the chimeflash range can be programmed to trigger one receiver.
- Works on the 868MHz frequency
- Range of up to 200m
- Easy to pair & completely Wireless



Sarabec Chimeflash Doorbell

| | | | |
|------------|---|--------|---|
| Order code | 150-003 | Colour | White |
| Price | £41.99 Ex. Vat £50.39 Inc. Vat | Size | Transmitter: 80x32x17mm Receiver: 145x105x52mm |



Chimeflash Bell Push

The Sarabec Chimeflash Doorbell Bell Push is a replacement or extra push that will pair to your Chimeflash Receiver unit. CR2032 Battery.

- Up to 200m range in open space
- Mounting plate with sticky pad
- Bright green indication light
- Battery expectancy: Up to 2 years



Chimeflash Telephone Transmitter

Get alerted when your telephone is ringing via extra loud sounds, ultra-bright flashing LED's and coloured indicators on your Chimeflash Receiver.



Chimeflash Bell Push

| | | | |
|------------|---|--------|-------------------|
| Order code | 150-021 | Colour | White |
| Price | £19.99 Ex. Vat £23.99 Inc. Vat | Size | 80x32x17mm |




Chimeflash Telephone Transmitter

| | | | |
|------------|---|--------|-------------------|
| Order code | 150-024 | Colour | White |
| Price | £27.65 Ex. Vat £33.18 Inc. Vat | Size | 80x32x17mm |



Chimeflash Sound Sensor

The new, easy-to-pair, Chimeflash Sound Sensor will attach to any device which emits a loud repeating sound such as a doorbell, door entry system, or alarm unit. It will send a wireless signal to the Chimeflash Receiver (sold separately with our Chimeflash Doorbell), alerting you via extra loud sounds, ultra-bright flashing LEDs, and coloured indicators.



Chimeflash Power Supply

The power supply for the Chimeflash Receiver System.

- 12V 1A output
- 1.5m Lead



Chimeflash Sound Sensor

| | | | |
|------------|---|--------|-------------------|
| Order code | 150-028 | Colour | White |
| Price | £27.65 Ex. Vat £33.18 Inc. Vat | Size | 80x32x17mm |



Chimeflash Power Supply

| | | | |
|------------|------------------------|--------|--------------|
| Order code | 150-049 | Colour | Black |
| Price | £17.99 Inc. Vat | | |



HONEYWELL SERIES 3 & 9 RANGE

The simple, versatile, and cost effective Series 3 range of door chimes from Honeywell are designed to be simple, contemporary, reliable, and hassle-free.

Including Honeywell ActivLink - the kits come pre-linked and ready to go straight out of the box.

Uncomplicated, screw-free installation means they are ready to go in minutes. Available in portable and plug in models.



FEATURES AND BENEFITS:

- Sleep/Mute mode silences the doorbell for either 3, 6, 9, or 12 hours or indefinitely
- Customisable halo light with 7 colour options
- Halo light & LED strobe indicators for visual alert
- USB charging capability allows you to charge your phone or tablet from the plug in modules.
- Adjustable maximum volume of 80-84 dB
- Choice of 4 or 6 crystal clear, digitally-enhanced chime sounds.
- Reliable wireless range of up to 150m*
- Wirefree; uses 868 MHz radio frequency to send signals between the bell push and the chime
- Up to 2 years battery life in portable units
- Lightweight and portable
- Easy-open battery compartments
- Secret Knock function - personalise melody & colours for friends/family
- Alarm accessories available to create a handy, simple entry alarm system.

**Max range between push and chime. Note: Building structures, e.g. metal and thick walls may affect transmission range.*





HONEYWELL PORTABLE DOORBELLS

Honeywell DC915N Halo Light Chime

- 7 colours and LED Strobe for visual alerts
- Wireless range of up to 200m
- 8 melodies
- Secret Knock function personalise melody & colours for friends/family
- Wall mounted OR portable chime
- Maximum volume of 90dB
- Extended battery life of up to 5 years



Honeywell DW915S Halo Light Chime Receiver

- Receiver only (no bellpush included)
- 8 melodies
- Battery powered or can be connected to an 8V/12V transformer
- 7 selectable colours and LED Strobe
- Link a second doorbell & expand the range up to 400m
- Adjustable volume with Sleep/Mute modes



| Honeywell DC915N Halo Light Chime | |
|---|--------------------------|
| Includes: Bellpush, Chime, CR2032 & C batteries | |
| Order code 212-920 | Colour White/Grey |
| Price £64.79 Inc. Vat | Size 126x126x42mm |
| | |

| Honeywell DW915S Halo Light Chime Receiver | |
|---|--------------------------|
| Additional Receiver. Includes: Chime, C Batteries | |
| Order code 212-923 | Colour White/Grey |
| Price £46.79 Inc. Vat | Size 126x126x43mm |
| | |

Honeywell DC315N Halo Light Chime

- Stylish white chime with Halo light
- White LED strobe indicators on top
- Choice of 6 melodies
- Up to 2 years battery life
- Easy installation - fixings included



| Honeywell DC315N Halo Light Chime | |
|---|--------------------------|
| Includes Chime, 4xAAA batteries, bell push & CR2032 | |
| Order code 212-320 | Colour White/Grey |
| Price £52.79 Inc. Vat | Size 110x70x43mm |
| | |

Honeywell DC315NBS Halo Plug-In Chime

- Stylish white plug-in chime
- White LED strobe indicators on top
- Choice of 6 melodies
- Night-light mode
- USB 2Amp charging port for mobile phones etc.



Honeywell DCP311 Bell Push Button

- Additional bell push for second door
- Weatherproof IP55 bell push with red LED confidence light
- Easy-open battery compartment



| Honeywell DC315NBS Halo Plug-In Chime | |
|---|--------------------------|
| Includes Chime & white bell push & CR2032 coin cell | |
| Order code 212-323 | Colour White/Grey |
| Price £58.19 Inc. Vat | Size 110x70x64 mm |
| | |

| Honeywell DCP311 Bell Push Button | |
|-----------------------------------|--------------------------------|
| Order code 212-331 | Colour White 70x30x16mm |
| Price £26.99 Inc vat | |
| | |

SIGNOLUX DOORBELL SYSTEMS

Alerting - Telling you what's happening



Signolux Tower Doorbell System

Introducing the new standard system in the Signolux range, designed for people with hearing problems to inform about domestic events. The 2-part system, made up of a Signolux Bellpush and Signolux Tower Receiver, comes supplied with rechargeable batteries and a mains power supply as standard.

- 16 x chime sounds
- Extra loud: 95dB+ on some sounds.
- Mains powered & rechargeable
- 2 x 3.5mm sockets for vibrating pad & other alerts
- 8 x easy-to-see visual indicators
- Bright flashing light with 360° visibility
- Compatible with all Signolux equipment



Signolux Tower Doorbell System

Includes: Flash & Chime Tower Receiver, Bellpush, Rechargeable Batteries & Charger.

| | | | | |
|------------|---|--------|-----------------|--|
| Order code | 230-750 | Colour | White | |
| Price | £72.49 Ex. Vat £86.99 Inc. Vat | Size | 143x80mm | |

Signolux Bell Push / Person Call Button

A vertical door bell push for the Signolux system which can be used at either the front or back door, or works well as a personal call button. The push features a red confidence light and can be wall-mounted using the supplied double-sided adhesive fixing strip.



Signolux Bell Push / Person Call Button

Bell push with CR2032 battery

| | | | | |
|------------|---|--------|-------------------|--|
| Order code | 230-710 | Colour | White/Grey | |
| Price | £16.15 Ex. Vat £19.38 Inc. Vat | Size | 80x40x16mm | |

Signolux Bell Push Landscape / Person Call Button

The Signolux wireless landscape bell push can be used as an additional door bell push for a Signolux visual alerting system or as a replacement.

- Weather resistant to IP55
- LED confidence light
- Name plate



Signolux Bell Push Landscape / Person Call Button

Landscape bell push with CR2032 battery and name card

| | | | | |
|------------|---|--------|-------------------|--|
| Order code | 230-711 | Colour | Grey/White | |
| Price | £16.15 Ex. Vat £19.38 Inc. Vat | Size | 80x40x16mm | |

Telephone (voice): 01642 247789
Fax: 01642 230827 Textphone (NGTS): 18001 07497 242 294
www.sarabec.com
Email: mail@sarabec.co.uk



AMPLICALL SENTINEL 1 VIDEO BABY MONITOR WITH LOUD ALARM & VIBRATING PAD



**£148.99
Ex. Vat**

Introducing our newest addition; the Amplicall Sentinel 1 Video Monitor includes a Video Monitor, a baby camera unit and a wired vibrating pad.

The vibration feature & bright flash act as an extra nudge for those with visual and hearing impairments.

When the baby monitor is triggered, the receiver will ring, flash, vibrate & display the Live feed onto your high-quality video receiver so you can sleep easy knowing your baby is too. Comes with a snake connector, which allows for very easy attachment of the camera unit to a baby cot.

- Camera with night mode
- Motion detector
- Temperature detector
- Noise detector
- Feeding time reminder
- LED indicators
- Connect 4 cameras to the same monitor
- Chat/zoom button



| Baby Monitor | | Colour | White |
|--------------|-------------------------------------|--------|-------|
| Order code | 252-120 | | |
| Price | £135.00 Ex. Vat £162.00 Inc. Vat | | |
| | | | |



ALARM CLOCKS

PORTABLE ALARM CLOCKS & TIMERS

TravelTim Portable Alarm Clock

Our most popular and very effective portable vibrating travel alarm clock.

- Digital alarm clock with backlit display
- 12/24 hr operation
- Strong vibration pattern
- Light and buzzer alarms
- Snooze control
- Pillow clip to secure clock to bedding
- 1xAAA clock battery and 2xAA vibrator batteries supplied
- Socket for connecting external vibrating pad (Order Code: 230-315)



TravelTim Portable Alarm Clock - Pink

| | | | | |
|------------|------------------------|--------|---------------------|---|
| Order code | 230-306 | Colour | Black/Pink |   |
| Price | £32.49 Ex. Vat | Size | 95 x 70 x 30 | |
| | £38.99 Inc. Vat | | | |

TravelTim Portable Alarm Clock - Green

| | | | | |
|------------|------------------------|--------|---------------------|---|
| Order code | 230-305 | Colour | Black/Green |   |
| Price | £32.49 Ex. Vat | Size | 95 x 70 x 30 | |
| | £38.99 Inc. Vat | | | |

iLuv Smartshaker 3

This vibrating bed shaker is a Bluetooth alarm clock that's ideal for heavy sleepers.

- 5 vibration strengths
- Multiple alarms
- Notifications & alerts
- Wake up in silence - no noise, just powerful vibrations

iLuv Smartshaker 3

| | | | | |
|------------|------------------------|--------|-------------------------|---|
| Order code | 230-339 | Colour | Black |    |
| Price | £49.99 Ex. Vat | Size | 43 x 135 x 175mm | |
| | £59.99 Inc. Vat | | | |



VIB Timer

A pocket-sized vibrating timer with a light pulse. Free-standing or clipped on, can be attached unobtrusively to your clothes – perfect for cooking, lab work, taking medication or when on the move.

- A compact timer with vibration, buzzer, and light signals
- Signal types/time units can be selected or switched on/off
- Timer function - counts down
- Stopwatch function



Vibrating Timer

| | | | |
|------------|------------------------|--------|--------------------------|
| Order code | 230-300 | Colour | |
| Price | £28.19 Inc. Vat | Size | Silver 70x45x20mm |



Sonic Shaker Portable Alarm Clock

The Sonic Shaker is a digital vibrating travel alarm clock. The clock is very easy to set up and use; with vibrate only, beep only, or vibrate with beep modes.

- 12-hour clock with PM indicator
- 4-minute snooze
- Clip & lead to secure in/under pillow
- Manual display LED light (when display is raised)
- Protective travel pouch for easy packing



Sonic Shaker Portable Alarm Clock

| | | | |
|------------|---|--------|-----------------|
| Order code | 205-100 | Colour | White |
| Price | £33.32 Ex. Vat £39.98 Inc. Vat | Size | 90Øx30mm |



Geemarc Wake 'n' Shake Voyager

This extra loud travel alarm clock with vibration is perfect for your travels. This compact, battery-operated alarm rings at a louder volume than a standard alarm, vibrates & transmits visual and audio alerts to wake even the deepest sleepers - effective for individuals with severe hearing impairment.

- Extra loud alarm (up to 75dB)
- Built-in vibrator with two levels of vibration
- 5x alarm settings
- Lamp and timer functions
- Easy-to-read display



Geemarc Wake 'n' Shake Voyager

| | | | |
|------------|---|--------|-------------------|
| Order code | 252-091 | Colour | White |
| Price | £28.99 Ex. Vat £34.79 Inc. Vat | Size | 95x90x30mm |



VC10 Vibrating Alarm Clock

A conventional analogue alarm clock with separate pad that has excellent vibration.

- Analogue alarm clock
- Illuminating clock face.
- Strong vibration pattern from slim, battery-powered pad
- Vibrate and/or beep alarm
- Snooze control
- Simple time and alarm controls
- 1AA clock battery and 3AA vibrator batteries supplied
- Can be used with additional flash/sound alerts



VC10 Vibrating Alarm Clock

| | | | |
|------------|---|--------|---|
| Order code | 230-310 | Colour | White |
| Price | £37.99 Ex. Vat £45.59 Inc. Vat | Size | Clock: 85x120x50mm Pad: 30x85x85mm |



MAINS-POWERED DIGITAL ALARM CLOCKS



Geemarc Wake 'n' Shake Curve

| | | | | |
|------------|------------------------|--------|---------------------------|--|
| Order code | 252-082 | Colour | Black | |
| Price | £44.99 Ex. Vat | Size | Clock: 178x94x50mm | |
| | £53.99 Inc. Vat | | Pad: 90Øx25mm | |

Geemarc Wake 'n' Shake Curve

The Wake 'n' Shake Curve is the very latest extra loud vibrating alarm clock from Geemarc. Benefiting from a new sleek curved design and a handy USB port for you to charge your smartphone or other devices.

- USB port to recharge mobiles and other electronic devices
- Extra loud (90dB) alarm
- Duration of alarm: 1 hour
- 4 alarm options
- Audible alarm with flashing LED light and vibrating shaker

Geemarc Wake 'n' Shake Star

A stylish, extra loud, extra bright alarm clock with vibrating pad that is certain to wake even the deepest of sleepers or those with hard of hearing. 5 alarm options, 2 vibration levels & a choice of 3 volume/tone levels.

- Extra loud alarm tone – 95dB
- White, extra bright flashing light (12 LEDs)
- Time display of 12/24 hours
- Lamp function for additional functionality
- Large white display with dimmer switch
- Can be used to charge smartphones and other electronics
- Back-up battery power (battery not included)



Geemarc Wake 'n' Shake Star

| | | | | |
|------------|------------------------|--------|-----------------|--|
| Order code | 252-084 | Colour | White | |
| Price | £45.83 Ex. Vat | Size | 90Øx30mm | |
| | £55.00 Inc. Vat | | | |



SB500 Sonic Bomb Alarm Clock

Do you struggle to get up in the morning? This Sonic Bomb alarm clock with a bed shaker is the perfect solution for heavy sleepers.

Stealth grey colour with red controls, night vision display, pulsating alert lights, and alarm ring indicator won't keep you in the dark. Complete 113db alarm with vibrating pad. Featured battery backup - maintains time and alarm setting in the event of power failure (battery not supplied).

SB500 Sonic Bomb Alarm Clock

| | | | | |
|------------|------------------------|--------|---------------------|--|
| Order code | 205-180 | Colour | Graphite/Red | |
| Price | £45.82 Ex. Vat | Size | 140x140x70mm | |
| | £54.98 Inc. Vat | | | |

- Adjustable volume & tone
- Snooze button/Alarm off button
- Dimmer switch
- AC power input
- Hour/minute button

MAINS POWERED ANALOGUE ALARM CLOCKS

Geemarc Wake 'n' Shake Vintage

This analogue, extra loud (80dB) alarm clock has a continuous 4-minute snooze, flashing ring light, and 5 x alarm settings. 5 types of alarm settings include: Off / flashing light display / vibrator & flashing light display / ringer & flashing light display / ringer & vibrator / flashing light display.

- Volume (High/Mid/Low/Off)
- Tones control (High/Mid/Low)



Geemarc Wake 'n' Shake Vintage

| | | | |
|------------|---|--------|---------------------|
| Order code | 252-092 | Colour | Black/White |
| Price | £49.99 Ex. Vat £59.99 Inc. Vat | Size | 140x140x70mm |



TIMEFLASH, VC10, & TRAVELTIM ACCESSORIES

Vibrating Pillow c/w Batteries

This Vibrating Pad can be connected to any of the Signolux receivers. This equipment uses its internal 3x AA batteries so does not draw power from the connected receiver, which is especially useful if using the system in a portable environment.



Vibrating Pillow c/w batteries

Comes with 1.5m lead.

| | | | |
|------------|---|--------|-------------------|
| Order code | 230-315 | Colour | Black |
| Price | £18.49 Ex. Vat £22.19 Inc. Vat | Size | 82x82x30mm |



Flash Module MF-1

Emits signals via powerful flashes of light, they make the whole room brighter and can be detected even with your eyes closed. Can be connected to any receiver from the Signolux Range.



Flash Module MF-1

| | | | |
|------------|---|--------|--------------------|
| Order code | 230-317 | Colour | Silver |
| Price | £56.17 Ex. Vat £67.40 Inc. Vat | Size | 65x93x134mm |



Sound Module MA-1

This module can be added to any existing Signolux receiver to allow you to hear the alerts as well as see and feel them. Generates a powerful alarm tone - continuously variable tone and volume settings.



Sound Module MA-1

| | | | |
|------------|---|--------|--------------------|
| Order code | 230-318 | Colour | Silver |
| Price | £49.81 Ex. Vat £59.77 Inc. Vat | Size | 134x93x65mm |



PHONES AND ACCESSORIES

CALLERS' VOICE TOO QUIET OR MUFFLED?

Do callers' voices seem very quiet or muffled over the telephone? In this catalogue, you'll find a selection of telephones; some that have been specially designed for people with hearing loss and others that have features that make them suitable for anyone with hearing impairments. Whether you're at home, in the office, or out & about, our selection of desktop (corded), portable (cordless), & mobile phones will suit your needs.

AMPLIFICATION AND TONE CONTROL

All telephones in this catalogue have variable amplification and most have a tone control. A tone control is extremely useful when trying to improve clarity, boost the frequencies that are being lost, and achieve clearer conversations. A mute function allows you to have a private conversation at your end without the person on the phone overhearing your conversation.

CORDED PHONES

All the landline corded telephones in this catalogue have an inductive coupler in the handset that can be used in conjunction with a hearing aid that has the "T" Loop setting or Loop Program. Some of the portable and mobile phones can also be used directly with hearing aids. If phones have this facility, they are listed as 'Hearing Aid Compatible' (HAC). Should interference be experienced, then a neckloop or ear hook can be used to eliminate the interference. Our Telephone Amplifiers for corded phones can be connected to almost all desktop phones to boost the volume and the tone. These amplifiers can be used at home or in the office.

PORTABLE OR DECT PHONES

Portable or DECT phones allow you the freedom to have a phone anywhere around the house or in the garden. Many have good amplification and tone controls. Some have an answer phone that you can access through the handset and benefit from the amplification. They can be used with hearing aids without interference.

MOBILE PHONES

Our mobile phones are easy to use and are accessible by all, with amplification and hearing aid compatibility as well as a strong vibration. We have mobile phones that are designed specifically for older users, which are simpler and don't have all the bells and whistles of some of the more complicated devices.





MOBILE PHONES

JUST SO YOU KNOW...

SIM CARDS

All our mobiles are unlocked (not linked to any specific network) and can be used with most SIM cards; either Contract or Pay As You Go. SIM cards are not supplied with our phones but are widely available to suit your locality. If you have an existing phone, you may be able to simply swap the SIM into your new phone and keep your old phone number.

BLUETOOTH®

Bluetooth® technology allows connections between different devices without the need for wires, and has a useful range of up to 10m.

HEARING AID COMPATIBILITY

Please remember that some hearing aids are more immune to interference than others. The amount and type of wireless signals transmitted by the phones can significantly vary, making it difficult to state which aids will work with which phones. However, we have found that the newer the hearing aid and the newer the phone, the less likely you are to have problems. You may find that you don't even need to use your "T" Loop setting or Loop Program - simply place the phone speaker as close to the hearing aid microphone as possible. A set of ratings exist to help users find mobile telephones that are suited for use with hearing aids.

M-RATINGS

Wireless devices rated for hearing aids switched to the microphone setting. The higher the number the better the rating. M4 is currently the highest.

T-RATINGS

Wireless devices rated for hearing aids switched to the T/Loop/telecoil setting. The higher the number the better the rating. T4 is currently the highest.

AMPLIFICATION

Manufacturers often quote a dB amplification range for their phones, this is **not** the absolute output of a handset, it is the volume control range of the phone from a level that is not quoted. It is just a guide to output.



DORO MOBILE PHONES

LOOKING FOR AN EASY-TO-USE MOBILE PHONE FOR YOURSELF OR FOR AN ELDERLY RELATIVE?

All our Doro phones are designed with seniors in mind and feature extra loud and clear sound, large separated keys, and a high visual contrast that make them easier to use than other phones.

Many Doro phones also come with a **'Response by Doro'** feature. This unique

service enables trusted relatives or friends to quickly and easily manage important settings on your phone, all via an app installed on theirs.

These services also keep your support network close. Just one press of your phone's assistance button sends an alarm with your location to a list of people you rely on, sounding a loud, distinct signal on their phones even if placed in 'Do not disturb' mode. More peace of mind for you, and for them.

Doro 8100 Smartphone

An easy-to-use smartphone that makes staying connected with family and friends both fun and reassuring. Conversations are clearer with Doro's optimized sound quality, and the intuitive interface simplifies everything you want to do.

- Intuitive menus
- Response by Doro
- Improved sound
- Triple camera
- Water resistant
- Face recognition
- Talk time: 23 hours
- Standby time: 450 hours

| Doro 8100 Smartphone | | Colour | Black |
|----------------------|---|--|-------------------|
| Order code | 253-343 | Size | 156x73x9mm |
| Price | £174.99 Ex. Vat £209.99 Inc. Vat |    | |



Doro 780x

Our **easiest-to-use mobile phone** that makes staying safe a whole lot simpler.

The water-resistant, easy-to-use Doro 780X has three large keys with name labels to put the people you rely on within easy reach without having to navigate menus or browse through a phonebook. It also has excellent safety functions such as a safety timer and an assistance button for sending an alarm with GPS position if you need help.

Numbers can be programmed remotely by a trusted relative via *Response by Doro* to automatically dial-up and send an SMS alarm with GPS location to preset numbers.

| Doro 780x | | Colour | White/Black |
|------------|---|---|--------------------|
| Order code | 253-102 | Size | 135x60x12mm |
| Price | £129.16 Ex. Vat £154.99 Inc. Vat |    | |



Doro 5860 white

The Doro 5860 is a great phone for easy everyday calling and texting. It features a wide landscape display, a camera for taking and sharing pictures, and an assistance button for alerting trusted contacts if you should ever need help.

| Doro 5860 white | |
|-----------------|---|
| Order code | 253-345 |
| Price | £66.66 Ex. Vat £79.99 Inc. Vat |
| Colour | White |
| Size | 66x13x128mm |



Doro 6820

The Doro 6820 is an excellent choice for seniors who want ease of use and great sound in a stylish clamshell design.

Calling, texting, and taking pictures are all made easier thanks to well-separated, high-contrast keys and a large 2.8" screen.

Other features include clearer conversations with HD voice, soft-touch coating for a more secure grip, and convenient shortcut keys to both camera and messages.

There is also an assistance button that can be used to alert up to five trusted friends or relatives if help is needed. And of course, it is compatible with hearing aid devices (HAC).

| Doro 6820 | |
|------------|---|
| Order code | 253-346 |
| Price | £83.32 Ex. Vat £99.98 Inc. Vat |
| Colour | White |
| Size | 107x56x20mm |



Doro 1380 Red

The Doro 1380 is a 2G phone which is great for people who prefer making calls, sending texts, and taking the occasional photos.

Dialling numbers and typing texts on the clear and spaced buttons is easy, while viewing messages or photos on the large, clear display is simple and enjoyable.

| Doro 1380 Red | |
|---------------|------------------------|
| Order code | 253-344 |
| Price | £34.99 Inc. Vat |
| Colour | Red |
| Size | 127x66x13mm |



PORTABLE DECT PHONES

Notes...

DECT

Digital Enhanced Cordless Telecommunications (Digital European Cordless Telecommunications), usually known by the acronym DECT, is a standard term primarily used for cordless phone systems.

GAP

A GAP (Generic Access Profile) compliant handset from one manufacturer should work, at a basic level, with a GAP-compliant base from another manufacturer, although it may be unable to access advanced features of the base station such as phone book synchronization or remote operation of an answering machine.



SWISSVOICE XTRA CORDLESS PHONES

SWISSVOICE XTRA ARE EASY-TO-USE, AMPLIFIED PHONES.

All Swissvoice phones featured here come with an answering machine and a call-block facility. Ideal for the hard of hearing and visually impaired individuals due to the loud volumes and large bright display. They also come with 4 direct memories function to easily call your favourite correspondents (with pictures on the corded base and A, B, C, D keys on the handset).

ADDITIONAL FEATURES:

- Hearing Aid Compatible (HAC)
- Vocalised handset
- Extra-loud ringer and volume boost
- Powerful ringers can be set up to 90dB
- Amplification of up to +40dB
- Large keypad with talking digits to help avoid dialling errors.
- Add up to 3 additional handsets
- Phonebook - up to 100 contacts

Swissvoice Xtra 2355 White
Cordless phone



Swissvoice Xtra 2355 Duo White
Cordless twin handset phone



Swissvoice Xtra 2355 White

| | | | |
|------------|---|--------|--|
| Order code | 255-621 | Colour | White |
| Price | £66.66 Ex. Vat £79.99 Inc. Vat | Size | Base: 170x118x60mm Handset: 185x54x25mm |



Swissvoice Xtra 2355 Duo White

| | | | |
|-------|--|--------|--|
| Code | 255-622 | Colour | White |
| Price | £91.66 Ex. Vat £109.99 Inc. Vat | Size | Base: 170x118x60mm Handset: 185x54x25mm |



Swissvoice Xtra 3355 Combo White Telephone

Duo phone with a corded base and a cordless DECT handset. Ideal to equip the whole house, this DECT combo combines the safety of a cordless phone with the flexibility of a cordless handset!

- High-quality, adjustable, hands-free function
- Ring indicator lights
- Ability to receive calls from stored numbers only

Swissvoice Xtra 3355 Combo White Telephone

| | | | |
|------------|--|--------|--|
| Order code | 255-627 | Colour | White |
| Price | £99.99 Ex. Vat £119.99 Inc. Vat | Size | Desktop: 217x200x82mm Cordless: 185x54x25 |



SWISSVOICE

Swissvoice 2355 Xtra Handset White

An extra handset for use or to be added to the Swissvoice range of Xtra cordless and combo phones in the 2355 & 3355 range.



Swissvoice 2355 Xtra Handset White

| | | | |
|------------|---|--------|-----------------------------|
| Order code | 255-629 | Colour | White |
| Price | £33.33 Ex. Vat £40.00 Inc. Vat | Size | Handset: 185x54x25mm |



DORO PHONEEASY DECT PHONES



PhoneEasy 100w DECT

An easy-to-use cordless phone with easy-to-read display, clear sound, and large clear buttons for easy operation.

- Amplification of up to +30dB
- Audio boost button
- Visual ring call indicator
- Tone control
- Hands-free speakerphone
- GAP compatible
- Up to 5 handsets per base

PhoneEasy 100w DECT

| | | | |
|------------|---|--------|--------------------|
| Order code | 253-069 | Colour | White |
| Price | £41.66 Ex. Vat £49.99 Inc. Vat | Size | 180x50x30mm |



PhoneEasy 100w DECT Duo

Two easy-to-use cordless phones with easy-to-read display, clear sound, and large clear buttons for easy operation. See above specs for '253-069: PhoneEasy 100w DECT'.

PhoneEasy 100w DECT Duo

| | | | |
|------------|---|--------|----------------------------------|
| Order code | 253-070 | Colour | White |
| Price | £58.32 Ex. Vat £69.98 Inc. Vat | Size | Handset: 180x 50x30mm |





GEEMARC AMPLIDECT 295 DECT PHONES

GEEMARC AMPLIDECT 295 TELEPHONES

If you find yourself straining to hear telephone conversations, then a handset offering adjustable amplified volume settings could help make life easier for you. The Geemarc AmpliDECT 295 Amplified Cordless Telephones with Answering Machine provides volume adjustment to help enhance your telephone calls.

The Geemarc Amplidect 295 Telephones are hearing aid compatible, and have:

- Amplification of up to +30dB
- Shaker function & adjustable ringer volume
- Built-in answering machine - 15 minutes of recording time and remote access to messages.
- Speakerphone
- Phonebook up to 50 numbers
- Redial function for last 10 outgoing calls
- P/T dialling mode
- R button for call waiting
- Hotkey to phonebook

Geemarc AmpliDECT 295

Helps enhance telephone use with volume adjustment and built-in answering machine.



Geemarc AmpliDECT 295

| | | | |
|------------|---|--------|--|
| Order code | 252-040 | Colour | White |
| Price | £58.32 Ex. Vat £69.98 Inc. Vat | Size | Handset: 172x55x25mm Base: 136x91x105mm |



Geemarc AmpliDECT 295 Twin

Twin handset version of the Geemarc Amplidect 295.



Geemarc AmpliDECT 295 Twin

| | | | |
|------------|--|--------|--|
| Order code | 252-041 | Colour | White |
| Price | £91.66 Ex. Vat £109.99 Inc. Vat | Size | Handset: 172x55x25mm Base: 136x91x105mm |



Geemarc AmpliDECT 295 Combi

A cleverly designed telephone system that combines the benefits of a photo ID phone with the flexibility of a cordless handset. Made up of a desktop corded telephone and a cordless handset that can be used anywhere around the house.

The desktop phone has big buttons and a large, easy-to-read display that is ideal for those with visual impairments, while the volume and tone is fully adjustable for people that are hard of hearing.

- 3 customisable one-touch photo buttons
- Large screen
- White backlight
- Alarm
- Possibility to pair up to 4 bases and 5 handsets
- Caller ID



Geemarc AmpliDECT 295 Combi

| | | | |
|------------|--|--------|--|
| Order code | 252-044 | Colour | White |
| Price | £91.66 Ex. Vat £109.99 Inc. Vat | Size | Base: 136x91x105mm Handset: 172x55x25mm |



Geemarc AmpliDECT 295 Extra Handset

You can pair up to 4 additional handsets to your AmpliDECT 295 and enjoy telephone conversations anywhere in your home; the range goes from 30 metres indoors to an incredible 80 metres outdoors!

Hear the difference with the Geemarc AmpliDECT 295 - a phone with more than just a ring! If you rely on the "T" Loop setting or Loop Program of your hearing aid to use your cordless telephone, then the Geemarc AmpliDECT 295 is the ideal solution for you! Comes with Caller ID.

Geemarc AmpliDECT 295 Extra Handset

| | | | |
|------------|---|--------|-------------------|
| Order code | 252-042 | Colour | White |
| Price | £33.32 Ex. Vat £39.98 Inc. Vat | Size | 172x55x25m |



Geemarc PhotoDECT 295 Extra Handset

Extra Handset for 295 systems.

Compatible with:

- AmpliDect295
- AmpliDect295 Twin
- AmpliDect Combi 295

Geemarc PhotoDECT 295 Extra Handset

| | | | |
|------------|---|--------|--------------------|
| Order code | 252-043 | Colour | White |
| Price | £37.49 Ex. Vat £44.99 Inc. Vat | Size | 172x55x25mm |



GEEMARC AMPLIDECT 595 DECT PHONES

GEEMARC AMPLIDECT 595 TELEPHONES

The Geemarc AmpliDECT 595 phones are ideal for anybody with hearing loss, dexterity or visual impairments.

This range offers large backlit keys for easy dialling, an adjustable receiving volume of up to +50dB, and a phonebook that can retain up to 50 names and numbers.

SOS FUNCTION BUTTON

For added security, you can save up to four contacts that can be alerted in times of emergency.

The cordless handset has a T coil, making it hearing aid-compatible for anyone with the "T" Loop setting or Loop Program activated on their hearing aid. Geemarc AmpliDECT 595 phones have:

- Answering machine & 30 mins recording time
- 10 one-touch memory buttons & pictures
- Amplification of up to +50dB
- Extra loud ringtones
- Voice announcement dialling
- Clock 12/24 hours
- Record voice message for memory photo buttons
- Visual ring indicator
- Hands-free speaker phone
- Automatic hang up when placed into base
- Redial function for last 10 outgoing calls
- Block up to 20 unwanted numbers
- Caller ID



Geemarc AmpliDECT 595

The Geemarc AmpliDECT 595 is ideal for anybody with hearing loss, dexterity, or visual impairments.

Geemarc AmpliDECT 595

| | | | |
|------------|--|--------|---|
| Order code | 252-045 | Colour | White |
| Price | £91.66 Ex. Vat £109.99 Inc. Vat | Size | Handset: 72x55x25mm Base: 136x105x91mm |



Geemarc AmpliDECT 595 Twin

The Geemarc AmpliDECT 595 TWIN is a dual handset version of the popular 595, ideal for anybody with hearing loss, dexterity, or visual impairments.

Geemarc AmpliDECT 595 Twin

| | | | |
|------------|---|--------|-------------------|
| Order code | 252-046 | Colour | White |
| Price | £133.32 Ex. Vat £159.98 Inc. Vat | Size | 72x55x25mm |



Geemarc AmpliDECT 595 Combi

Ideal for anybody with hearing loss, dexterity, or visual impairments.

Geemarc AmpliDECT 595 Combi

Order code **252-048**
Price **£133.32 Ex. Vat**
£159.98 In. Vat

Colour **White**



Geemarc AmpliDECT 595 / 595 U.L.E. Extra Handset

An extra/additional handset to compliment the Geemarc 595 range of phones.

Geemarc AmpliDECT 595 / 595 U.L.E. Extra Handset

Order code **252-049**
Price **£58.32 Ex. Vat**
£69.98 Inc. Vat

Colour **White**
Size **72x55
x25mm**



Geemarc AmpliDECT 595 U.L.E.

Ergonomic and easy-to-use cordless handset, with a big button keypad and backlight.

With the Doorbell 595 U.L.E. accessory, that has the ability to forward the communication onto your mobile or another corded telephone, you will be able to speak with the person at your door when you are not home.

- VIP memory allows incoming calls only for the numbers saved in the VIP phonebook (up to 50 numbers)
- Handset LCD display
- White backlight

Geemarc AmpliDECT 595 U.L.E.

Order code **252-053**
Price **£91.66 Ex. Vat**
£109.99 Inc. Vat

Colour **White**
Size **Handset: 185x55x25mm**
Base: 136x91x105mm



CORDED DESKTOP PHONES



Telephones (Desktop/Corded)

AMPLICOMMS POWERTEL 196



The Amplicomms PowerTel 196 is an extra-loud corded, amplified telephone, designed with hearing impaired and/or elderly people in mind who often need an extra-loud phone to hear conversations clearly. A lot of phones focus solely on volume, quite often at the expense of sound quality. This phone has a volume increase of +60dB but also includes in-built technology that makes sure the reproduced sound is crisp and clear.

- Two dedicated speed-dial keys (M1 and M2)
- Four-line extra-large point matrix display screen with blue backlight for easier viewing



- Compatible with hearing aids (HAC)
- Amplification up to +60dB
- Eight handset volume levels can be easily adjusted
- Audio boost key for users with severe hearing loss
- Extra-loud ringer (90dB)
- Extra-bright ring indicator
- Large keys with large easy-to-see numbers
- Built-in tone equaliser with 3 settings
- Visual alert for incoming calls
- Caller ID presentation (CLIP) with capacity to store the last 30 callers
- Telephone book for storing up to 200 numbers
- Redial function for the 10 most recent numbers
- Suitable for both domestic and business use

Amplicomms PowerTel 196

Order code **255-202**
 Price **£66.58 Ex. Vat**
£79.90 Inc. Vat

Colour **Grey/Silver**
 Size **195x220x108mm**



SWISSVOICE XTRA 1110 WHITE



This easy-to-use corded phone is made for people needing maximum comfort. Pin the 6 photo buttons with your relatives' pictures to instantly call them. In addition, the phone is fitted with 2 direct memory keys and is Hearing Aid Compatible.

- 2-position cursor on audio boost (On/Off)
- 1 ringtone
- 2 melody levels + Off
- Max melody volume - 85dB at 50 cm
- Voice EQ on headphone: Bass, Treble, Natural
- Possible wall-mounting
- Marking type FV/DC

The handset offers 4 levels of volume, whilst the audio-boost key offers an extra loud earpiece volume of +30dB. To avoid missing any calls, the ringer can be set to +85dB (+ off position). The phone has a large backlit incoming call indicator, and features a hands-free function (8 levels) with a backlit key for the pleasure of sharing conversations.

This is a standalone phone, not compatible with Xtra cordless handsets



- Hearing aid compatible
- Amplification up to +30dB
- Hands-free function with backlit key
- 6 direct memories (M1, M2)
- Keyboard with large keys with Braille writing
- Large incoming call indicator light
- Last number redial key
- Chat
- Adjustable volume on earphone

Swissvoice Xtra 1110 White

Order code **255-600**
 Price **£33.33 Ex. Vat**
£40.00 Inc. Vat

Colour **White**
 Size **217x200x82mm**



DORO CORDED PHONES

Magna 4000

Amongst the loudest amplified telephones with excellent layout and features. One push of the boost button instantly raises the handset volume. Very easy to use thanks to large keys, direct memories, and a display for Caller ID. Includes phonebook with 50 names and numbers and other useful features.

- Amplification of up to +60dB
- Loud ringer with two settings
- Adjustable tone control
- Call indicator light
- 2 Number quick dial
- Wall mount option
- Caller ID with 50 number recall
- Hands-free speakerphone
- Inductive coupler for hearing aid users



PhoneEasy 311c

The Doro 311c is an easy to use large button telephone with a comfortable handset and an adjustable volume control.

- Amplification of up to +11dB
- Inductive coupler for hearing aid users with the "T" loop setting or program
- Visual ring indicator
- 3 ringer melodies
- 3 direct memories
- 10 indirect memories
- Last number redial



Magna 4000

Order code **253-026**
Price **£70.82 Ex. Vat**
£84.98 Inc. Vat

Colour **Black**
Size **200x220x85**



PhoneEasy 311c

Order code **253-011**
Price **£29.99 Inc. Vat**

Colour **White**
Size **185x220x90**



PhoneEasy 312cs

The Doro 312cs is an easy-to-use telephone with a large good contrast display, hands-free speakerphone and visual ring indicator. It has a 30-number phonebook & visual display for Caller ID, dialled numbers, and phonebook.

- Amplification of up to +18dB
- Inductive coupler for hearing aid users with the "T" loop setting or program
- Adjustable ring volume to 85dB with 9 ringtones
- 3 direct memories
- Last number redial with 15 memories



PhoneEasy 331 Photo

The PhoneEasy 331 is an easy-to-use, large button telephone with clear, well-spaced & high contrast keys, and 3 photo memories.

- Amplification of up to +18dB
- Inductive coupler for hearing aid users with the "T" loop setting or program
- Visual ring indicator
- 3 direct photo memories
- High contrast keys with low reflective finish
- Last number redial



PhoneEasy 312cs

Order code **253-012**
Price **£44.99 Inc. Vat**

Colour **White**
Size **185x220x90mm**



PhoneEasy 331 Photo

Order code **253-031**
Price **£32.99 Inc. Vat**

Colour **White**
Size **185x230x100mm**



GEEMARC CORDED PHONES



Geemarc CL100 Big Button Telephone

| | | | |
|------------|---|--------|---------------------|
| Order code | 252-171 | Colour | White |
| Price | £39.99 Ex. Vat £47.99 Inc. Vat | Size | 195x210x70mm |



Geemarc CL100 Big Button Telephone

The CL100 is a perennial favourite amongst amplified corded telephones. With a large button keypad, ringer volume control, and the ability to amplify incoming speech. Fully hearing aid compatible.

- Amplification of up to +30dB
- Amplifies outgoing speech for those with a quiet voice 4dB
- Tone control for incoming speech with good boost for high frequency loss +/- 10dB
- Inductive coupler for hearing aid users with the "T" loop setting or program
- Call indicator light
- 3 emergency & 9 memory number buttons
- 4 AAA Batteries (not supplied) for memory backup
- Wall-mountable

Geemarc CL1100

The CL1100 is a standard desktop amplified corded telephone with good amplification, including the ability to amplify incoming and outgoing speech. Fully hearing aid compatible.

- Amplification of up to +30dB
- Outgoing speech volume control +4 dB
- Variable tone control +/-10dB
- Inductive coupler for hearing aid users with the "T" loop setting or program
- Visual ring indicator
- Adjustable ringer level, tone & melody
- 12 one-touch direct memories
- Last number redial & Mute function



Geemarc CL1100

| | | | |
|------------|---|--------|---------------------|
| Order code | 252-173W | Colour | White |
| Price | £39.16 Ex. Vat £46.99 Inc. Vat | Size | 170x220x50mm |



Geemarc CL595 Amplified Telephone

| | | | |
|------------|---|--------|---------------------|
| Order code | 252-219 | Colour | White |
| Price | £82.29 Ex. Vat £98.75 Inc. Vat | Size | 220x200x70mm |



Geemarc CL595 Amplified Telephone

The Geemarc CL595 is an extra loud, corded phone with large photo buttons, an SOS function, and an answering machine. It has raised tactile keys, perfect for users with dexterity issues. Fully hearing aid compatible.

- Amplification of up to +50dB
- 4 one-touch memories (with talking caller ID), 3 with pictures
- Speakerphone
- Redial function for last 10 outgoing calls
- Phonebook can store up to 50 names and numbers (4 of which have talking caller ID)
- Mute/Pause function

Geemarc AmpliPower 40

One of the telephones with a high level of amplification, the AmpliPOWER 40 is a similar design to the CL100. It also has a bright flashing call indicator and greater choice of ringer volume and sounds. Designed to help those with considerable hearing loss and visual impairments.

- Amplification of up to +40dB
- Volume boost button
- Tone control for incoming speech with good boost for high frequency loss +/- 10dB
- Inductive coupler for hearing aid users with the "T" loop setting or program
- Hands-free speakerphone
- Ringer volume and tone choice controls
- 12 Number Memory
- Socket (2.5mm) for headset, neckloop, or ear hook use
- Socket for optional vibrating pad



Geemarc AmpliPower 40

| | | | |
|------------|---|--------|---------------------|
| Order code | 252-194 | Colour | White |
| Price | £62.49 Ex. Vat £74.99 Inc. Vat | Size | 195x210x70mm |



Geemarc AmpliPower 50

| | | | |
|------------|---|--------|---------------------|
| Order code | 252-195 | Colour | White |
| Price | £83.32 Ex. Vat £99.98 Inc. Vat | Size | 195x210x70mm |



Geemarc AmpliPower 50

Amongst telephones with the highest amplification the AmpliPOWER 50 is similar to the AmpliPOWER 40, but it is the loudest and clearest amplified telephone available, with a maximum 60dB amplification. It is designed to help those with a severe to profound hearing loss. Compatible with all hearing aids that have a telecoil facility.

- Amplification of up to +60dB
- Volume boost button
- Receiving tone control up to 10dB
- Inductive coupler for hearing aid with the "T" loop setting or program
- Hands-free speakerphone
- Bright call indicator light
- Ringer volume and tone choice controls
- Socket (2.5mm) for headset, neckloop, or ear hook
- 12 Number Memory
- Headset and shaker outputs

Geemarc Serenities Telephone

The Geemarc Serenities big button amplified phone is ideal for people with hearing or visual impairment. With amplification up to 30dB and extra bright visual ring indicator. SOS pendant included - it allows the user to call for help with one press of the pendant button, which automatically dials up to 5 emergency contact numbers and plays a personalised voice announcement upon pickup.

- Hearing Aid Compatible (for hearing aids set to the T-setting)
- Amplification of up to +30dB
- Adjustable receiving tone control (4 levels)
- Outgoing speech amplification (-4, 0, & +4dB)
- 5 one-direct memories and 10 priority memories
- Hands-free/speakerphone
- Redial key
- Mute function
- Phonebook (30 numbers)
- Waterproof (IP67) emergency transmitter included (range of 20 to 30m)
- Wristband auto dials up to 5 emergency numbers
- Confirmation tone call



Geemarc Serenities Telephone

| | | | |
|------------|---|--------|--------------------|
| Order code | 252-218 | Colour | White |
| Price | £116.66 Ex. Vat £139.99 Inc. Vat | Size | 20x180x50mm |





PHONE ACCESSORIES



PhonePlus Telephone Handset

A battery-operated telephone amplifier for most desktop, corded handset telephones, which works by increasing the volume of the caller's voice through their handset, while boosting high/low tones if needed. Simply disconnect the handset from the phone base and connect the PhonePlus. Features LED Indicator and Tone Control.

Note: Not suitable for use with cordless, mobile or telephones with buttons in the handset, & not ideal for already highly amplified phones.

- Volume control range of up to 40dB with boost button
- Volume control range of +/- 10dB
- Socket for optional headset, neckloop, or ear hook
- Battery powered; 2 alkaline AA supplied, approx life 2 years
- 20cm connection lead for phone comes supplied

PhonePlus Telephone Handset

| | | | | |
|------------|---|--------|--------------------|--|
| Order code | 140-051 | Colour | Silver | |
| Price | £29.99 Ex. Vat £35.99 Inc. Vat | Size | 102x70x35mm | |



iBell2 Silver Mobile Phone Alert

iBell2 informs you about missed calls and shows every ring of your mobile phone with blue running lights. The iBell2 is triggered by the vibration of the mobile phone. The vibration alerting for incoming calls and messages must be activated within the phone settings.

Should you have a Humanteknik alerting system, you can also include your mobile phone within the system. The iBell2 can be connected via a Signalux Acoustic Universal Transmitter. Suitable for all mobile phones and pads with a vibrating option.

iBell2 Silver Mobile Phone Alert

| | | | | |
|------------|---|--------|------------------|--|
| Order code | 230-374 | Colour | Silver | |
| Price | £39.90 Ex. Vat £47.88 Inc. Vat | Size | 120x145mm | |



Geemarc iLoop+ Neckloop with Microphone

Powered neckloop to connect to any Smartphone, iPhone, MP3, Kindle etc. Connects to any audio device with a 3.5mm headset socket & 2 x AAA (1.2V rechargeable NiMH batteries with up to 2 weeks standby time).


- Switch between music and phone calls at the touch of a button
- Built-in microphone
- Voice Activation and auto shut off when not in use

Geemarc iLoop+ Neckloop with Microphone

| | | | | |
|------------|---|--------|---------------------|--|
| Order code | 252-187 | Colour | Black/Silver | |
| Price | £54.16 Ex. Vat £64.99 Inc. Vat | Size | 75x35x15mm | |

TINNITUS

Explore products to help relieve symptoms of **tinnitus** and sleep soundly with relaxing background noises.



To explore our bluetooth sleep headphones & sound players such as the Snoozebands, you can turn to [page 16](#) of this catalogue.



Tinnitus

Tinnitus

Wellcare Tinnitus Sound Relaxer

This product can help relieve the symptoms of tinnitus through relaxation or distraction. It acts as a stress-reliever, helping you achieve a restful night's sleep.

- Choice of 7 sounds: heartbeat, sea waves, running stream, birdsong, woodland glade, falling raindrops or simple white noise
- Volume control and separate on/off switch
- 30-minute sleep timer
- Battery powered - 4 x AA batteries



Wellcare Tinnitus Sound Relaxer

| | | | |
|------------|---|--------|--------------------|
| Order code | 258-005 | Colour | White/Grey |
| Price | £20.99 Ex. Vat £25.19 Inc. Vat | Size | 140x110x130 |



Sound Oasis Sound Therapy System S-680-02

| | | | |
|------------|---|--------|---------------------|
| Order code | 258-105 | Colour | Black |
| Price | £62.49 Ex. Vat £74.99 Inc. Vat | Size | 120x87x143mm |



Sound Oasis Sound Therapy System S-680-02

Benefit from specially developed sounds for tinnitus relief as well as the ability to customise the sound list to your own preference.

- World-leading sound quality – sleep & nature sounds play at a higher quality than a CD to lull you to sleep
- Seamless Sleep-Enhance experience – Sounds gradually slow, allowing to you gently slip into a restful sleep
- Custom sound cards – select a mix of up to 24 of Sound Oasis' expertly crafted sleep and relaxation sounds
- Improved clock and alarm setting features
- Refreshed & new sleek look
- USB powered – misplaced your adapter? Power this white noise machine with a standard USB adapter

Sound Oasis Travel Sleep Sound Therapy (S-850W)

This unique, compact sound therapy device combines high-performance digital sounds with a full-featured travel/alarm clock. Help manage tinnitus by blocking out annoying hotel noises & creating a soothing sleep environment while travelling or at home.

Drift off to sleep with soothing nature sounds and wake up to one of 4 alarm types: sound, buzzer, voice memo, or the exclusive jet lag reduction sounds. On/off functions include: A unique "gradual off" volume reduction, a rotary on/off volume control, & a one-touch off/resume button.



- 18 digitally recorded sounds
- 19 authentic, high fidelity sounds
- Voice memo - record your own alarm or memo
- 2 x time zone jumbo backlit alarm clock & calendar
- 4 x position timers: choose 30, 60, or 90 minutes
- Headphone/speaker jack
- Travel pouch & time zone chart included
- Powered by 4 "AA" batteries or AC adaptor (included)

Sound Oasis Travel Sleep Sound Therapy

| | | | |
|------------|---|--------|--------------------|
| Order code | 258-107 | Colour | Silver |
| Price | £66.55 Ex. Vat £79.86 Inc. Vat | Size | 137x38x76mm |



Telephone (voice): 01642 247789
Fax: 01642 230827 Textphone (NGTS): 18001 07497 242 294

www.sarabec.com
Email: mail@sarabec.co.uk



Sound Oasis Bluetooth Speaker BST-100

Sleep soundly with any one of the 10 built-in sounds to help relieve tinnitus and to provide a calming background noise for relaxation and contemplation. With Bluetooth capabilities, this speaker gives the option to connect to a mobile phone or computer wirelessly and stream additional sounds, music, or even an audio book.

- Perfect for improving sleep, blocking noise, managing tinnitus, and calming babies.
- Pre-loaded microSD card with 10 superior quality 'Sound Oasis' sounds.
- Your selected sound plays seamlessly for 8 hours.
- Exceptional big bass sound, but super compact size is perfect for use at home or while travelling.
- Stream new sounds or music from any Bluetooth-enabled device for exceptional versatility. Alternatively you can connect your audio device (MP3 player, CD player, mobile phone etc) via an auxiliary cable (not supplied).
- Insert a new microSD card with 'Sound Oasis' sounds for unlimited sound library expansion.
- 30'/10m Bluetooth wireless range.
- Built-in rechargeable lithium battery allows for 5 - 8 hours of cordless use on one charge.



Sound Oasis Bluetooth Speaker BST-100

Order code **258-109**
 Price **£44.80 Ex. Vat**
£53.76 Inc. Vat

Colour **Silver/Grey**
 Size **63Øx66mm**



Sound Oasis Stereo Pillow Speakers SP101

Order code **258-120** Colour **White**
 Price **£20.39 Inc. Vat** Size **126x63x20mm**



Sound Oasis Stereo Pillow Speakers SP101

A pair of small stereo pillow speakers with handy in-line volume control from Sound Oasis.

The low profile pillow speakers simply slip under your pillow to allow you to listen to any radio, CD, or cassette player, etc., that has a 3.5mm stereo speaker/headphone socket. They are wired together so you don't lose one, but you are still able arrange them comfortably under your pillow.

Use with Sound Oasis Sound Relaxers, CD players, iPods, iPhones, iPads, most MP3 players, and TVs.

- Two ultra-thin speakers for comfortable use
- In-line volume control for effortless volume control.
- High quality stereo sound
- Speaker output: 4 ohm, 2-watt speakers
- Speaker Jack: 1/8" (3.5 mm) plug fits most electronic devices including Sound Oasis Sound Therapy Systems, iPods/MP3 players, CD players, TVs, etc

COMMERCIAL AUDIO

IMPROVING COMMUNICATION FOR THE HEARING IMPAIRED

We can work with you to improve communication for hearing impaired people throughout various, non-domestic environments.

We offer systems that meet the continually varying needs of customers to be compliant with current legislation in order to provide quality audio services for providers, clients and users.

APPLICATIONS

- Reception areas
- Conference and meeting rooms
- Courtrooms
- Lecture theatres
- Churches
- Multilingual simultaneous translation
- Tour guides
- Auditoriums
- Multimedia rooms and theatres

BESPOKE SERVICE

We can specify and design a bespoke system that fits your type of building and its requirements. Through our Loop System Design package, we can deliver a full supply, install, and provide an after-care package. We can work with you to improve communication for hearing impaired people throughout various community settings.

BENEFITS OF LARGE AREA LISTENING

The main advantage for all systems is that the sound is delivered to the listeners in a pure, undistorted form that is not influenced by the distance between the listeners and the sound source, or by any annoying background noise in the room.

LETS GET IT RIGHT

In order to find the right listening system for your needs, it is important to compare the requirements of the venue and its architectural features to ensure that the properties of the chosen transmission system is not in conflict with its surroundings. Frequent site visits and

assessments are an important part of any pre-installation for a new or upgrading venue.


The key to all these transmission systems is their input from an audio source, such as the microphone; one or many, or a public address system. If the quality of the input is not at its best, then the output will not be satisfactory.

HOW THE SYSTEMS WORK

An induction loop system consists of a loop of insulated wire (the loop) that runs around the required listening area, an amplifier and microphone, or other input. The amplifier transmits the required audio signals as an alternating current into the loop. A magnetic field of audio sound is then created inside the loop, allowing listeners with the "T" loop setting or loop program on their hearing aids, to move freely around the room and receive the audio signals directly into their hearing aids without unwanted sounds.

INSTALLATIONS

Please contact us if you have a project in mind and would like further information or help with a system, system design, installation, commissioning, and training.





LA90 PORTABLE INDUCTION LOOP SYSTEMS

The LA90 Portable Loop Systems provide a convenient solution for hearing aid users to listen more easily to conversation in noisy situations via the "T" Loop setting or Loop Program of their hearing aid. Used in reception areas, meeting rooms, interviews, checkouts, etc. for a 1:1 conversation when in noisy areas or where discretion is required.

A built-in microphone within the unit picks the sound up of the person speaking and transmits a loop signal to the hearing aid user. The speaker and hearing aid user must be within approximately 1m(3ft) of the unit. The hearing aid user will see a confidence light on the front of the unit to show them the system is operating, while the speaker will see an indication of their voice being transmitted. The Loop System transmits a signal both ways so both individuals can receive the signal equally well.

In a group setting, the system can be placed flat on a table, by the separate microphone, in front of the speaker.



LA90 Desktop Loop System

This system can be used with an external desktop microphone to enhance a person's speech. It also sends signals directly to the hearing aids & gives more discretion.

| LA90 Desktop Loop System | |
|---|---|
| Includes: LA90 Loop System, Desktop Microphone, Handset | |
| Order code | 137-102 |
| Price | £190.90 Ex. Vat £229.08 Inc. Vat |
| Colour | Silver |
| Size | 200x185x70 |
| | |

LA90 Portable Induction Loop

A Portable loop amplifier with internal microphone, external microphone, and headphone socket (for users without the "T" Loop setting or Loop Program), volume control, replaceable NiMH battery pack, and mains power supply. The additional microphone picks up voices from both users and transmits clearly with less background noise. Battery life of approx. 6 hours from a 4-hour charge.

- Facility to fix amplifier to counter top or use Kensington lock
- Ear loop window stickers to advertise loop service



| LA90 Portable Induction Loop | |
|------------------------------|---|
| Order code | 137-100 |
| Price | £116.65 Ex. Vat £139.98 Inc. Vat |
| Colour | Silver + choice of 2 panel colours |
| Size | 200x185x70 |
| | |



LA90 Handset Kit

This adds a telephone handset to the standard LA90, enabling hearing aid and non-hearing aid users alike to benefit from easy to access amplification.

| LA90 Handset Kit | |
|------------------|---|
| Order code | 137-106 |
| Price | £32.49 Ex. Vat £38.99 Inc. Vat |
| Colour | Graphite |
| Size | 200x185x70 |
| | |

PROLOOP SYSTEMS



The ProLoop Amplifier family is designed for professional use, covering any room from small to very large, giving reliable dependable use where its most required.

ProLoop systems incorporate Automatic Gain Control (AGC), which guarantees consistent field strength and reproduces stable sound with a special emphasis on clear speech, even under acoustically demanding conditions. Corresponding test equipment can be used to evaluate the sound quality via the integrated monitor out jack. With the ProLoop systems, you'll receive high availability, operating reliability, and quality assurance.

The **ProLoop NX3** and **NX7** systems are primarily designed for continuous use. The integrated switch mode power supply and the Class-D amplifier design enable high efficiency and low waste heat. Among other things, this results in high reliability and operational safety. The function of the induction loop and the amplifier is continuously monitored. A functional error is indicated via displays on the unit and at the error output. Both systems incorporate Automatic Gain Control (AGC) and Metal Loss Correction (MLC), which can be used to correct the frequency response when affected by metal structures present in the building.



ProLoop C (170m2)

The ProLoop is a switchable, 100% short-circuit-proof amplifier, with two symmetrical XLR input jacks, and an extremely stable (4 A RMS) output performance.

- Automatic fuse reset & extra-sturdy jacks
- 1x RCA input jack & 1x Line output jack
- Dual action AGC for excellent speech recognition
- Convenient monitoring of the magnetic field (can be done using headphones or speakers)
- Treble control to compensate for treble losses due to reinforcement



ProLoop C (170m2)



| | | | |
|---|-------------------------|--------|---------------------|
| Order code | 230-200 | Colour | Black |
| Price | £335.88 Inc. Vat | Size | 64x295x205mm |
|   | | | |

ProLoop NX3

- Designed for small-sized rooms (up to 56m²).
- The balanced input 1 can be switched between microphone & line sensitivity, & provides 48V phantom power for condenser microphones
- The unbalanced Input 2 handles line level signals or the audio signal of an electret microphone
- Input 3 is designed for connecting media players via a 3.5mm stereo jack. Stereo signals are summed to mono.



ProLoop NX3

| | | | |
|---|-------------------------|--------|--------------------|
| Order code | 230-209 | Colour | Black |
| Price | £250.80 Inc. Vat | Size | 32x167x95mm |
|   | | | |

ProLoop NX7

- 7A RMS Output Current
- Link connections allow several PROLOOP NX Amplifiers to be combined into one unit
- The output signal can be monitored via a headphone output
- A screw-fixed cover plate prevents settings from changing
- The power switch on the back prevents the unit from being switched off unintentionally

**ProLoop NX7**

Order code **230-213**
Price **£598.80 Inc. Vat**

Colour **Black**
Size **43 x430x290mm**



ACCESSORIES FOR COMMERCIAL AUDIO

**Crescendo 60/2 System**

A system which can be used as a personal listener as well as to listen to, or monitor existing loop systems. Supplied with Under-chin Stethoscope Headset, Plug-in Microphone, External Microphone (1.2m lead), 3m Extension Lead, and 2 x AA batteries.

**Crescendo 60/2 system**

Order code **109-102**
Price **£79.99 Ex. Vat**
£95.99 Inc. Vat

Colour **White/Grey**
Size **95x62x25mm**



Sarabec

**ProLoop Field Strength Meter Plus**

The Proloop Field Strength Meter (FSM) Plus measures the magnetic field strength of induction loop systems according to IEC 60118-4:2006. The device is designed for professional use and delivers reliable RMS values (RMS-measuring range calibrated to 400mA/m = 0dB) on the output level, response frequency, AGC function, distortion and background noise in the tested induction loop. In addition, it enables the user to make an acoustic evaluation of the sound via earphones. It offers adjustment possibilities for test ranges of background noise and field strength, and multi-coloured LEDs for excellent legibility.

- LED-scale with 3dB resolution
- Handy, lightweight, rugged chassis
- Switch-able A weighted filter
- Uses two 1.5V AA batteries & has 3.5mm stereo headphone jack with volume control

ProLoop Field Strength Meter Plus

Order code **230-211**
Price **£130.79 Inc. Vat**

Colour **Grey**
Size **83x36x125mm**



Ordering Information

ORDERING INFORMATION

You can view and order all of our products online at www.sarabec.com 24 hours a day 7 days a week. We use a secure ordering system to ensure the safety of your details.

If you wish to place your order by post or fax, simply tear out the order form overleaf and complete it with the details (product code, price, and quantity) of the products you require.

- If you are eligible for VAT relief, then, for each product required, you only pay the Excluding VAT price - You will need to complete the VAT Relief Certificate on the back of the order form.
- If a VAT exclusive price is not quoted for the product then the price including VAT must be paid.
- If paying by VISA, Mastercard, Maestro or AMEX Please complete the credit/debit card payment details.
- Post or fax your order to us.
- If paying by cheque, make it payable to 'Sarabec Ltd'. Enclose your payment with this order and post to us. Please note that cheques are cleared before goods are sent.

VAT

Throughout this catalogue you will find products that are priced Inclusive (Inc. VAT) and Exclusive (Ex. VAT) of VAT. VAT relief is not available on postage and packing.

If, in the description of each product we have detailed prices for both 'Including' and 'Excluding' VAT, that product is eligible for VAT relief. If an exclusive price is not shown, VAT relief is not available.

You are entitled to VAT relief if you are handicapped by a hearing disability. To benefit, simply complete the VAT Relief Certificate included in this catalogue. If you are unable to make this declaration, you must include the VAT in your order.

PAYMENT

You can pay by using a cheque or postal order payable to SARABEC Ltd or by using a credit card or debit card. If you wish to use a credit or debit card you may place your order and pay online via secure payment services or direct by telephone, fax or textphone.

TRIAL PERIOD

All products carry a one month trial period if purchased within the United Kingdom. You may return your products to us for any reason within one month of receipt and obtain a refund for your purchase provided that:

- You contact us first to let us know why you are returning the product
- The products are returned in good condition and in their original packaging; we may charge if products have been misused.
- The products have been used in accordance with the product instructions.

The cost of returning any product is your responsibility. We cannot accept responsibility for loss or damage to goods until they arrive safely at our offices. We advise you to return goods adequately packed and covered by insurance where appropriate.

DELIVERY

Products are normally available from stock and are despatched within 48 hours of receipt of your order. Delivery to UK mainland addresses normally takes up to **3 working days** from receipt of your order by us. Delivery is made through Royal Mail or independent carriers. A signature may be required on delivery of the product. Our delivery charges are for all UK mainland addresses, except the Channel Islands. For addresses outside this area please ask for delivery costs before ordering.

FURTHER INFORMATION

Some products in this catalogue may not be available for the whole life of the catalogue.

All personal information entered on the Order or VAT Relief forms will be treated as confidential as detailed in our Privacy Policy. You can ask for a copy of the Privacy Policy or obtain one from our website.

Our normal hours of business are: 8.15am to 4.45pm - Monday to Thursday, 8.15am to 4.00pm Friday.

**Sarabec Limited, 15 High Force Road,
Middlesbrough, TS2 1RH, United Kingdom**

Telephone: (voice) 01642 247789

Fax: 01642 230827

Textphone (NGTS): 18001 07497 242 294

Website: www.sarabec.com

E-mail: mail@sarabec.co.uk

PRICES

Whilst we do our best to keep any price rises to a minimum, we reserve the right to increase our prices at any time. Should we do so, we will contact you on receipt of your order with details of the increase giving you the opportunity to cancel your order if you wish.

Should your payment not clear, then we may, in addition to any other rights we may have:

- Suspend delivery of the product ordered or any part of them until we receive cleared funds.
- Legally recover those products already delivered to you from you or your premises at any reasonable time and retain ownership.
- Pass on to you any charges we incur in obtaining clearance of your payment.

You will not own the products until your payment for them has cleared.

VAT

VAT is calculated at the prevailing rate at date of invoice. Should the prevailing rate of VAT change from that used in our publications, calculations will be based on new rates and charged to you accordingly.

DELIVERY

Delivery to UK addresses normally takes 3 working days from the date of despatch of your order by us. It is made through independent carriers. The independent carrier may require a signature from you on receipt of your order.

If you do not receive your order within 28 days of placing your order or should your order arrive damaged, please contact us as soon as reasonably practicable. Please retain the original packaging.

Additional postage charges may be payable for deliveries to non-UK addresses; please contact us for further details.

GUARANTEE

All our products are covered by our guarantee. You will invalidate the guarantee if:

- You use the products contrary to instructions supplied with them
- You cause any damage to the products through

misuse or negligent use of them.

- If any product becomes damaged outside normal usage.
- For any consumables such as batteries.

If you wish to rely on this guarantee, we will, at our choice, either repair or replace the products or, if appropriate, any parts of them. If we are unable to carry out any of the above through no fault of our own we will refund the price you paid for the products.

In all circumstances, except for faulty goods, this guarantee shall not exceed the price you paid for the products in question. No other remedies are available.

OUR PRODUCTS

Full descriptions of our products are given in our current catalogue or online.

The products are only suitable for the intended uses described.

As a part of our continuing programme of product development, we reserve the right to improve and update product designs.

CANCELLATION

You may cancel your order anytime from when you place the order up to 10 working days after the date you received the goods but, to do so, you must notify us. We will respond confirming your cancellation and the date your notification was received. The date confirmed is the date the cancellation came into effect.

UNFORESEEABLE PROBLEMS

We are not responsible for delays outside our control (for example industrial dispute, act of God, inability to obtain raw materials) if this prevents or hinders our contract with you.

OUR AGREEMENT

An agreement is made between us once we have received and accepted your order.

Our agreement and any related matters are governed by English law. Your statutory rights are not affected. This agreement is a summary of our full Terms and Conditions of Sale and Supply; you can ask for a copy or obtain one from our website.

June 2018

ORDER FORM

Please complete the form below with the details of the products you wish to order.

If you are eligible for VAT relief, then, for each product required you only pay the price excluding VAT. If a VAT exclusive price is not quoted for the product then the price including VAT must be paid.

| Order Code | Description | Price Excluding VAT | Price Including VAT | Quantity | Total Price £ |
|---|-------------|---------------------|---------------------|----------|---------------|
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| Postage & packaging for UK mainland addresses ONLY: ■ Standard Carriage: £5.99 (Inc. Vat) Please enquire for delivery to areas outside this area ■ 24-Hour Carriage: £8.99 (Inc. Vat) | | | | | |
| Grand Total payable to 'Sarabec Ltd' (cheques may take up to 7 days to clear) | | | | | |

DELIVERY NAME & ADDRESS

If claiming VAT relief please complete the certificate overleaf.

| | |
|----------------|---------------|
| Name: | |
| Address: | |
| | |
| | |
| Postcode: | Telephone no: |
| E-mail address | |

If you have entered an email address above, please can we send your invoice to that address? YES / NO
Would you like to receive details of our products from us in future? (Please tick) By post? By email?

CREDIT/DEBIT CARD PAYMENT

I wish to pay by (please tick):    

Please debit my card number

| | | | | | | | | | | | | | | | | | | | |
|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
| | | | | | | | | | | | | | | | | | | | |
|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|

| | | |
|--|-------------------|------------|
| Card Valid from: | Card Expiry date: | Issue No.: |
| CVV Security Code: <i>The security code consists of the last three/four digits found on the signature strip.</i> | | |
| Card holder's name as shown on card (PLEASE PRINT): | | |
| Card holder's signature: | Date: | |
| Address: | | |
| | | |
| Postcode: | Telephone no: | |

To help us improve our service please can you tell us where you obtained this catalogue

VAT RELIEF CERTIFICATE

Please note there are penalties for making false declarations

If you are in any doubt as to whether you are eligible to receive goods or services zero-rated for VAT you should consult notice 701/7 reliefs for disabled people (viewable at www.hmrc.gov.uk) or contact the National Advice Service on 0845 010 9000 before signing the declaration.

| |
|------------|
| Full Name: |
| Address: |
| |
| |
| |
| Postcode: |

I declare that:

- I have a hearing disability
- I am receiving from Sarabec Ltd, 15 High Force Road, Middlesbrough, TS2 1RH the items overleaf which are being supplied to me for my domestic or personal use.
- I claim relief from VAT.

| |
|------------|
| Signature: |
| Date: |



SARABEC LIMITED

15 High Force Road, Middlesbrough TS2 1RH

Telephone (voice) 01642 247789 | Fax 01642 230827 | Textphone (NGTS): 18001 07497 242 294

Website: www.sarabec.com | E-mail: mail@sarabec.co.uk



Accredited Contractor
CHAS.co.uk



SARABEC LIMITED, 15 High Force Road, Middlesbrough, TS2 1RH, UK
TEL: 01642 247789 | **TEXTPHONE:** 18001 07497 242 294 | **FAX:** 01642 230827
EMAIL: mail@sarabec.co.uk | **WEBSITE:** www.sarabec.com